

INFINITY

COLLECTION

INFINITE POSITION SEATING

The Ultimate in Relaxation and Support

The Infinity Collection from Pride® Power Lift Recliners features infinite position seating and pocketed coils designed to deliver comfort, style and performance. Available in three different fabric selections, the Infinity Collection models deliver maximum relaxation and will match any décor.



Easy-to-use
hand control with
USB charging port



Scan here for more
information on the
Infinity Collection.



FOLLOW US ON SOCIAL MEDIA

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

FEATURES & BENEFITS

- Pocketed coils deliver individualized support and relaxation
- Arm and headrest cover standard.
- Premium furniture-grade laminate/hardwood frames for added strength and durability
- Easy-to-use hand control with USB charger (not available with optional heat and massage)
- Integrated battery backup ensures chair is operational when power goes out (batteries not included)
- Energy-efficient power box equipped with self-tested electronics and 9' cord
- Standard lumbar pillow
- Large, dual pockets for extra storage
- Quiet and smooth lift system

OPTIONS*

- Deluxe heat and massage
- Residential cleaning kit
- Head pillow
- Protective pad

* Additional charges apply for all options. Extra lead time may apply for some options.



MODEL	LC-525iS	LC-525iM	LC-525iL
POSITION	Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'3" & Below	Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'4" - 6'0"	Infinite Position Zero Gravity Trendelenburg Chaise Lounger 5'9" - 6'2"
SIZE	Small	Medium	Large
WEIGHT CAPACITY	375 lbs.	375 lbs.	375 lbs.
SEAT-TO-FLOOR	19"	20"	21"
SEAT DEPTH	19"	20"	22"
SEAT WIDTH	19"	20"	22"
TOP OF BACK TO SEAT	26"	28"	30"
BACK STYLE	Blown Fiber Split Back	Blown Fiber Split Back	Blown Fiber Split Back
FULLY PADDED CHAISE	Yes	Yes	Yes
DELUXE HEAT AND MASSAGE	Optional	Optional	Optional
SEATING TYPE	Pocketed Coil Springs	Pocketed Coil Springs	Pocketed Coil Springs
DISTANCE FROM WALL	20"	20"	21"
OVERALL WIDTH	34"	34.5"	37"
CHAIR WEIGHT	130 lbs.	130 lbs.	130 lbs.
WARRANTY	7-Year Pro-rated	7-Year Pro-rated	7-Year Pro-rated

DISCONTINUED

NEW FABRICS
COMING SOON!

Durasoft® Standard

Lexis Sta-Kleen®

UltraFabrics® Ultraleather®

Crypton® Aria

Beamsville, ON / 888-570-1113 / pridemobility.ca

Note: Measurements may vary +/- 1". The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.