

# Jazzy<sup>®</sup>

## SELECT 6




The Jazzy Select<sup>®</sup> 6 incorporates sporty design and performance features for a truly unique ride experience. Pride's patented Active-Trac<sup>®</sup> Suspension combined with six wheels on the ground and an innovative articulating beam allow for superior performance indoors and out, with added stability over uneven terrain and ramp transitions.

### ACTIVE-TRAC<sup>®</sup> SUSPENSION



**JAZZY ARMOR<sup>®</sup> AND  
COLOUR-THROUGH SHROUD**



**COMFORT HIGH-BACK SEAT  
WITH REPLACEABLE BACK  
AND SEAT COVERS**

**NORTH AMERICA'S  
FAVORITE POWER CHAIR**

Learn more online.



[www.pridemobility.ca](http://www.pridemobility.ca)

# Specifications

<b>Weight Capacity</b>	300 lbs.
<b>Maximum Speed</b>	Up to 4 mph
<b>Ground Clearance</b>	1.75"
<b>Turning Radius</b>	21.875"
<b>Overall Length</b>	35.375"
<b>Overall Width</b>	23.5"
<b>Drive Wheels</b>	10"
<b>Front Wheels</b>	6" solid casters
<b>Rear Wheels</b>	6" solid casters
<b>Suspension</b>	Active-Trac® (ATX)
<b>Drivetrain</b>	Two-motor, in-line, mid-wheel drive
<b>Braking System</b>	Regenerative and electro-mechanical
<b>Standard Electronics</b>	40 amp, PG GC 3 controller
<b>Battery Charger</b>	Off-board, 3.5 amp
<b>Per-Charge Range</b>	Up to 15 miles
<b>Battery Requirements</b>	(2) 12 volt, U-1
<b>Battery Weight</b>	25 lbs. each
<b>Base Weight</b>	100 lbs.
<b>Standard Seat Weight</b>	43 lbs. (comfort, high-back)
<b>Standard Seat</b>	16" x 16-18", 18" x 18-20", 20" x 18-20"
<b>Max. Seat Dimensions</b>	20" x 20"
<b>Warranty</b>	5-year limited on frame; 13-month limited on batteries; 13-month limited on electronics; 13-month limited on drive motors

Colours:



Jazzy Red



Jazzy Blue

# Features and Benefits

- Jazzy® Armor® includes shroud and controller guards to protect against daily wear and tear
- Depth adjustable, limited recline comfort high back seat available
- Universal, fixed, in-line joystick mount
- Flip-up, angle-adjustable foot platform
- Black, flat-free, non-marking tires
- 60" lap belt

## Seating Options\*

- Depth adjustable, limited recline comfort high back seat with solid seat pan, available sizes: 16"x16"-18", 18"x18"-20", 20"x18"-20"
- Deluxe contour, limited recline high back available sizes: 18"x16"
- Deluxe contour, limited recline, high back with solid seat pan, available sizes: 16"x16"-20"x20"
- Synergy® Seat
- PES STF range of 22.5"-26.5"

## Options

- Bell shaped foot platform
- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away footrests
- Swing-away joystick
- VR2 6-key (must upgrade with PES option)

## Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Lap belt (50", 70", 80", 90")
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger

\*Additional charges may apply.