

Model	Jazzy Passport
Part Number	FJPT1700
Weight Capacity	250 lbs.
Maximum Speed ¹	Up to 3.6 mph
Ground Clearance ³	3" at motor
Turning Radius ³	24.25"
Overall Length ³	40"
Overall Width ³	23.5"
Seat-to-Ground Height ³	19"
Seat-to-Foot Platform Height ³	14.5" (no adjustment)
Casters ³	8" solid
Drive Tires ³	11.25" perforated
Suspension	Limited
Range Per Charge ^{1,2}	Up to 7.9 miles (18 AH @ 250 lbs.) Up to 8.6 miles (18 AH @ 220 lbs.) Up to 4.6 miles (10.4 AH travel battery @ 250 lbs.) Up to 5 miles (10.4 AH travel battery @ 220 lbs.)
Total Weight Without Batteries ³	59.8 lbs. (unit does not disassemble)
Heaviest Piece When Disassembled ³	N/A (unit does not disassemble)
Standard Seating ³	Material: Black, Foam Dimensions: 17.5" x 16.5"
Drivetrain	Two-motor, rear-wheel drive
Dual Braking System	Regenerative and electromechanical
Standard Electronics	50 amp, PG VSI
Battery Requirements ³	(1) 25.2V 18 AH Lithium ion (1) 25.2V 10.4 AH Lithium ion
Battery Weight ⁴	7.6 lbs. (18 AH) 5.2 lbs. 10.4 AH Lithium ion (airline-compliant travel battery)
Battery Charger	2-amp, off-board
Waranty	Frame: 5-year limited Drivetrain: 24-month Electronics: 24-month Battery: 13-month

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Includes standard seat weight.

Jazzy[®] PASSPORT



NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power chair only. Overall measurements will vary based on seating and accessory selections.



Note: Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Jazzy[®] Passport and its lithium ion battery are permitted.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]

Beamsville, Ontario • 888-570-1113

pridemobility.ca |    

CANJAZZYPASSPORT_100121

Jazzy[®] PASSPORT

Simply Fold and Travel

Airline
Compliant
Travel Battery
Available



**NORTH AMERICA'S
FAVORITE POWER CHAIR**

pridemobility.ca |    

MAKES TRAVELING A BREEZE



Folds in just a few simple steps
for effortless transport

Experience the beauty of lightweight travel with the **Jazzy® Passport!** This compact power chair folds in just a few simple steps for effortless transport. With convenient storage the Jazzy Passport makes traveling a breeze.

- A compact and lightweight design
- Back seat pocket for extra storage while on the go
- Includes convenient, under seat mesh storage bag, mesh cup holder and 40" lap belt
- Large foot platform to accommodate the user's comfort
- Specialty foam seat design for extreme comfort
- Choice of 18 AH battery or 10.4 AH* travel battery
- Flat-free, perforated tires

Colours



Black

Accessories

- XLR USB charger



Optional

- Stealth Simplicity G Cushion, 18"x18"



Folded Size
31" x 16" x 23.5"

*Airline-compliant.