

END-TO-END SCOOTER SOLUTIONS



| | ZERO TURN SERIES | | FULL-SIZE LUXURY | | | TRAVEL MOBILITY | | | | | MID-SIZE LUXURY | VALUE SERIES | | HEAVY-DUTY | | | | |
|--|--|---|--|--|--|--|--|---|-----------------------------|---|-----------------|--|--|--|--|---|--|-----------------------------|
| | Zero Turn 10 ¹⁸ | Zero Turn 8 | Wrangler™ | Raptor ¹⁸ | Pursuit® XL ¹⁸ | Go-Go® iCRUISE | Go-Go® Endurance | Go-Go® LX with CTS Suspension ¹⁰ | Go-Go® Sport ¹⁰ | Go-Go® Folding Scooter ¹⁰ | Go-Go® Ultra X | Victory® LX with CTS Suspension ¹⁸ | Victory® 10.2 ¹⁰ | Revo® 2.0 ¹⁰ | Victory® Twin ¹⁸ | Maxima | Pursuit® ES ¹⁸ | |
| FRONT TIRES | 9" SOLID | 2" x 7" SOLID | 14.5" TUBELESS PNEUMATIC | 3" x 16" TUBELESS PNEUMATIC | 13" PNEUMATIC | 6.5" SOLID | 2" X 8" SOLID | 2" x 7" SOLID | 2" x 7" SOLID | 1.6" x 7" SOLID | | 3.5" x 10" SOLID (PNEUMATIC OPTIONAL) | 3.5" x 10" SOLID | 3" x 8" SOLID (3W) 3" x 9" SOLID (4W) (PNEUMATIC OPTIONAL) | DUAL 3" x 9" SOLID | 3W: 10.5"x3.6" SOLID 4W: 10" x 3.6" SOLID | 4" x 13" PNEUMATIC | |
| REAR TIRES | 10.75" SOLID | 2.5" x 9" SOLID | 14.5" TUBELESS PNEUMATIC | 3" x 16" TUBELESS PNEUMATIC | 13" PNEUMATIC | 6.75" SOLID | 2.5" X 8" SOLID | 2.5" x 8" SOLID | 2.5" x 8" SOLID | 2" x 8" SOLID | | 3.5" x 10" SOLID (PNEUMATIC OPTIONAL) | 3.5" x 10" SOLID | 3" x 9" SOLID (PNEUMATIC OPTIONAL) | 3.5" x 10" SOLID | 10.75" x 3.6" SOLID | 4" x 13" PNEUMATIC | |
| MAXIMUM SPEED ² | UP TO 7 MPH | UP TO 4 MPH | UP TO 11.4 MPH (200 LBS.) UP TO 11 MPH (350 LBS.) | UP TO 14 MPH | UP TO 7.5 MPH | UP TO 3.7 MPH | UP TO 4.2 MPH | UP TO 4.4 MPH | UP TO 4.4 MPH | UP TO 4 MPH | | UP TO 5.6 MPH | UP TO 5.2 MPH | UP TO 5 MPH | UP TO 5.25 MPH | UP TO 5.8 MPH | UP TO 6 MPH HS: UP TO 9.5 MPH | |
| GROUND CLEARANCE ^{4,13,17} | 2" AT MOTOR | 1.4" AT MOTOR | 5.5" AT MOTOR | 3.5" AT MID-FRAME | 5.25" AT MID-FRAME | 1.3" AT MOTOR | 3.25" AT SIDE DECK | 1.5" AT MOTOR | 1.5" AT MOTOR | 2.3" | | 2.4" AT TRANSAXLE | 2.4" AT TRANSAXLE 2.1" AT ANTI-TIP 3.0" AT SIDE DECK | 3W: 3.25" 4W: 3.5" AT MID-FRAME | 3" AT MID-FRAME | 3W: 2.5" AT MOTOR 4W: 2.25" AT MOTOR | 3.5" | |
| 3-W TURNING RADIUS ⁴ | N/A | N/A | N/A | 60" | N/A | 33.5" | N/A | N/A | N/A | N/A | | N/A | 42" | 40" | N/A | 49" | N/A | |
| 4-W TURNING RADIUS ⁴ | 43" | 38.25" | 94.75" | N/A | 82.5" | N/A | 48.25" | 44.75" | 44.75" | 32" | | 60" | 52.25" | 52" | 42" | 53" | 68.5" | |
| 3-W LENGTH/WIDTH ^{4,12} | N/A | N/A | N/A | 59"/27.13" | N/A | 38"/19.25" | N/A | N/A | N/A | N/A | | N/A | 45"/22.25" | 46"/20.75" | N/A | 47.5"/27" | N/A | |
| 4-W LENGTH/WIDTH ^{4,12} | 48"/24.25" | 45.75"/20" | 61.25"/32.5" | N/A | 56"/27.25" | N/A | 42"/20.5" | 41"/20.5" | 41"/20.5" | 37"/19" | | 47"/22.25" | 47.25"/22.25" | 47"/20.75" | 46"/26" | 47.5"/27" | 53.5"/25" | |
| SEAT-TO-GROUND HEIGHT RANGE | 22.25"-24.25" | 21.5"-23.5" | 30"-31.5" | 28" | 26.25"-27.75" | 18"-20" | 21.25"-22.75" | 20.75"-22.25" | 20.75"-22.25" | 20.5" | | 23.5"-25.5" | 22"-24" | 21"-23.5" | 23"-24.5" | 22.76"-24.76" (3W) 22.1" - 24.1" (4W) 22.32"-28.56" (with power elevating seat) | 24.75"-26.25" | |
| SEAT TYPE (STANDARD) | COMPACT, FOLDABLE (18" X 17") | COMPACT FOLDABLE, MOLDED PLASTIC (17"X17") | LIMITED RECLINING HIGHBACK PILLOW TOP (18"X16.5") | HIGH BACK (20" X 18") (100" LIMITED RECLINE) | HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18") | COMPACT (14.75" X 13.5") | COMPACT, FOLDABLE (17" X 17" AND 20" X 18") | COMPACT, FOLDABLE (17"X17") | COMPACT, FOLDABLE (17"X17") | COMPACT, FOLDABLE (16"x13") | | HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18" X 18") | FOLDABLE (18"X17") | COMPACT, FOLDABLE (20"X18") | FOLDABLE, MOLDED PLASTIC (18"X17") | DEPTH-ADJUSTABLE LIMITED RECLINE HIGH-BACK (22" x 18"-20") | HIGH-BACK RECLINER WITH SEAT SLIDE, SWIVEL, ADJUSTABLE ARMS AND HEIGHT-ADJUSTABLE HEADREST (18.5" X 18.5") | |
| (3/4-WHEEL) WEIGHT ^{3,4,14} | N/A / 185.8 LBS. | 120 LBS. | 274 LBS. | 208 LBS. | 259.9 LBS. | 42 LBS. | 93.5 LBS. | 85 LBS. | 85 LBS. | 51 LBS. | | N/A/156.25 LBS. | 130.25 LBS./139 LBS. | 120 LBS./132 LBS. | N/A/150 LBS. | 183.5 LBS. (3W) 196.5 LBS. (4W) | N/A/161 LBS. | |
| SEAT-TO-DECK HEIGHT RANGE | 15.5"-17.5" | 17"-18.5" | 18.75"-20.25" | 19" | 21"-22.5" | 8.75"-10.75" | 16.25"-17.75" | 16.25"-17.25" | 16.25"-17.25" | 13" | | 16.5"-18.5" | 15.25" - 17.25" | 15"-18" | 17.5"-19" | 3W: 16.25"-18.25" 4W: 17.13" - 19.13" 17.36"-23.61" (w/power elevating seat) | 17.75"-18.75" | |
| BATTERY WEIGHT ¹¹ | 30.25 LBS. EA | 13.5 LBS. EA. | 48 LBS. EA. (75 AH) 60 LBS. EA. (100 AH) | 53 LBS. EA. | 70 LBS. EA. (G31) 60 LBS. EA. (G27) 53 LBS. EA (G24) | 4.5 LBS. | 8 LBS. EA. (8 AH) 12 LBS. EA. (16 AH) | 13.3 LBS. EA. ⁸ | 13.3 LBS. EA. ⁸ | 9 LBS. EA. (SEALED LEAD ACID) / 5 LBS. EACH (LITHIUM-ION) | | 24 LBS. EA./ 32 LBS. EA. | 24 LBS. EA./ 32 LBS. EA. | 25 LBS. EA. | 24 LBS. EA./ 32 LBS. EA. | 36.5 LBS. EA. | 53 LBS. EA. | |
| BATTERY REQUIREMENTS ⁷ | (2) 40AH (INCLUDED) | (2) 18AH (INCLUDED) | (2) 75AH OR (2) 100 AH (OPTIONAL) | (2) 75AH (INCLUDED) | (2) GROUP 31 ¹⁵ , (110AH) (2) GROUP 27 (100AH) OR (2) GROUP 24 (75 AH) | (1) 36 VOLT, 5.2 AH LITHIUM ION (INCLUDED) | 8 AH LI-ION PACK OR 16 AH LI-ION PACK | (2) 18AH (INCLUDED) | (2) 18AH (INCLUDED) | (2) 12AH SEALED LEAD ACID (STANDARD) OR (1) 11.5 AH LITHIUM ION ¹⁶ | | (2) U-1 ¹⁶ OR (2) 40AH ¹⁶ | (2) U-1 (INCLUDED) OR (2) 40AH ¹⁶ | (2) U-1 (INCLUDED) | (2) U-1 (INCLUDED) OR (2) 40AH ¹⁶ | (2) GEL, 58AH, NF-22 | (2) GROUP 24 ¹⁶ | |
| PER CHARGE RANGE (MAX.) ^{2,9} | UP TO 24 MILES (200 LBS.) UP TO 18 MILES (400 LBS.) | UP TO 11.6 MILES (325 LBS.) UP TO 13.4 MILES (200 LBS.) | UP TO 18.6 MILES AT 350 LBS. (75AH); UP TO 24.2 MILES AT 350 LBS. (100AH); UP TO 22.6 MILES AT 200 LBS. (75AH); UP TO 29.4 MILES AT 200 LBS. (100AH) | 19 MILES (UP TO 400 LBS. AND FULL SPEED) 31 MILES (UP TO 220 LBS. AND AT LOW SPEED) | G31: UP TO 31.83 MILES G27: UP TO 29.51 MILES G24: UP TO 22.13 MILES | UP TO 10.3 MILES (250 LBS.) UP TO 12.7 MILES (200 LBS.) | UP TO 7.2 MILES (8 AH) UP TO 14.4 MILES (16 AH) | UP TO 14.6 MILES | UP TO 14.6 MILES | 4W: UP TO 13 MILES (LITHIUM-ION) 4W: UP TO 9.3 MILES (SEALED LEAD ACID) | | UP TO 16.16 MILES | 3W/4W: UP TO 11 MILES (U-1) 3W: UP TO 13 MILES (40AH) 4W: UP TO 13.5 MILES (40 AH) | 3W: UP TO 17.4 MILES (200 LBS.); 3W: UP TO 13.2 MILES (400 LBS.); 4W: UP TO 17.8 MILES (200 LBS.); 3W: UP TO 12.4 MILES (400 LBS.) | UP TO 11.1 MILES (WITH U-1 BATTERIES) | 3W: UP TO 19.8 MILES 4W: UP TO 19.5 MILES | UP TO 33.6 MILES HS: UP TO 22 MILES | |
| BATTERY CHARGER | OFF-BOARD, 3.5A | OFF-BOARD, 2A | OFF-BOARD, 8A | OFF-BOARD, 8A | OFF-BOARD, 8A | OFF-BOARD, 36 VOLT LITHIUM ION | OFF-BOARD, 3.5A LITHIUM-ION | OFF-BOARD, 2A | OFF-BOARD, 2A | OFF-BOARD, 2A (SEALED LEAD ACID) OFF-BOARD, 2.5A (LITHIUM-ION) | | OFF-BOARD, 5A | OFF-BOARD, 3.5A (U-1) OFF-BOARD, 5A (40 AH) | OFF-BOARD, 3.5A | OFF-BOARD, 3A (U-1) OFF-BOARD, 5A (40 AH) | OFF-BOARD, 5A | OFF-BOARD, 8A | |
| SUSPENSION | FRONT & REAR | FRONT & REAR | FRONT & REAR | FRONT & REAR | FRONT & REAR | N/A | FRONT & REAR | FRONT & REAR | FRONT & REAR | N/A | | FRONT & REAR | N/A | FRONT & REAR | N/A | N/A | FRONT & REAR | |
| MAXIMUM WEIGHT CAPACITY | 400 LBS. | 325 LBS. | 350 LBS. | 400 LBS. | 400 LBS. | 250 LBS. | 325 LBS. | 300 LBS. | 300 LBS. | 250 LBS. | | 400 LBS. | 400 LBS. | 400 LBS. | 350 LBS. | 500 LBS. | 350 LBS. | |
| STANDARD BODY COLOURS | CANDY APPLE RED, OCEAN BLUE, PEARL WHITE, BLACK | FIRE OPAL, GREEN MOONSTONE, ONYX BLACK, PEWTER, PEARL WHITE, SAPPHIRE BLUE ¹ | DESERT SAND, STONE GREY | BLACK (WITH COLOR INSERT OPTIONS) | CANDY APPLE RED, SILVER | LIME, MANGO, ARCTIC ICE, BLACKBERRY, RASPBERRY | GARNET RED (MATTE), EURO GRAY | RED, BLUE ¹ | RED, BLUE ¹ | WHITE | | CANDY APPLE RED, VIPER BLUE | CANDY APPLE RED, OCEAN BLUE | GREY STREET, TRUE BLUE | CANDY APPLE RED, VIPER BLUE | CANDY APPLE RED, VIPER BLUE | CANDY APPLE RED, VIPER BLUE | CANDY APPLE RED, VIPER BLUE |
| WEIGHT OF HEAVIEST PIECE ⁴ | 75 LBS. (FRONT SECTION) | 48 LBS. (REAR SECTION) | 219 LBS. (UNIT DOES NOT DISASSEMBLE) | 313.6 LBS. (UNIT DOES NOT DISASSEMBLE. WEIGHT INCLUDES SEAT, BASE AND BATTERIES) | 205.5 LBS. (MAIN FRAME) | 32 LBS. (SCOOTER BASE) | 40.75 LBS. (FRONT SECTION) | 37 LBS. (FRONT SECTION) | 37 LBS. (FRONT SECTION) | 51 LBS. (UNIT DOES NOT DISASSEMBLE - UNIT WEIGHT WITHOUT BATTERY PACK) | | 77.6 LBS. (FRONT SECTION) | 56.5 LBS. (3-WHEEL) (FRONT SECTION) 64.5 LBS. (4-WHEEL) (FRONT SECTION) | 48 LBS. (3-WHEEL) (REAR SECTION) 57 LBS. (4-WHEEL) (FRONT SECTION) | 55 LBS. (FRONT SECTION) | 76 LBS. (3-WHEEL) (REAR SECTION) 76 LBS. (4-WHEEL) (REAR SECTION) | 118 LBS. (MAIN FRAME) | |

FOOTNOTES

- All shroud panels included with choice of colour
- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- Battery weight not included.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of + or - 3%.
- 12 AH battery pack weighs 18 lbs.
- Durable solid plastic.
- Deep cycle (AGM or gel cell type recommended).
- 18 AH battery pack weighs 28 lbs.
- Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- Includes pathway illumination light.
- Battery weight may vary based on manufacturer.
- Width of base (not including seat).
- Does not account for anti-tip clearance.
- Includes standard seat weight.
- Option available through most providers. Not available through Pride Mobility[®].
- Optional.
- At midframe.
- Light package.

The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed & range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%.

ZERO TURN SERIES



ZERO TURN 10
4-wheel scooter

- Ocean Blue (Matte)
- Candy Apple Red
- Pearl White
- Black (Matte)

FEATURING
iTURN™
TECHNOLOGY



ZERO TURN 8
4-wheel scooter

- Sapphire Blue
- Pewter
- Pearl White
- Black
- Green Moonstone
- Fire Opal

Interchangeable
Shrouds

FEATURING
iTURN™
TECHNOLOGY

TRAVEL MOBILITY



GO-GO® iCRUISE
3-wheel scooter

- Lime
- Arctic Ice
- Mango
- Blackberry
- Raspberry

GO-GO® ENDURANCE
4-wheel scooter

Coming Soon



- Garnet Red (matte)
- Euro Gray



GO-GO® LX WITH CTS SUSPENSION
4-wheel scooter

OPTIONAL

- Red
- Blue
- Orion Yellow
- Onyx Black
- Pearl Pink
- Sunburst Orange

Interchangeable
Shroud Inserts

GO-GO® SPORT
4 & 3-wheel scooter

- Red
- Blue



GO-GO® FOLDING SCOOTER
4-wheel scooter

- White

Coming Soon



GO-GO® ULTRA X
4 & 3-wheel scooter

- Black

FULL-SIZE LUXURY

WRANGLER™
4-wheel scooter

- Desert Sand
- Stone Gray



RAPTOR
3-wheel scooter

- Black
- True Blue
- Candy Apple Red
- Green Machine
- Lemon Crush
- Pink Topaz
- Evolution Orange
- Matte Black

Optional
Interchangeable
Shroud Inserts



PURSUIT® XL
3-wheel scooter

- Candy Apple Red
- Viper Blue



MID-SIZE LUXURY



VICTORY® LX WITH CTS SUSPENSION
4-wheel scooter

- Candy Apple Red
- Viper Blue

VALUE SERIES



VICTORY® 10.2
4 & 3-wheel scooter

- Candy Apple Red
- Ocean Blue (Matte)



VICTORY® TWIN
4-wheel scooter

- Candy Apple Red
- Ocean Blue (Matte)



REVO® 2.0
4-wheel scooter

- True Blue
- Grey Street

HEAVY DUTY



MAXIMA®
4 & 3-wheel scooter

- Garnet Red (matte)
- Ocean Blue

ALSO AVAILABLE
IN HIGH-SPEED MODEL



PURSUIT® ES
4-wheel scooter

- Candy Apple Red
- Viper Blue