

## MODEL

Jazzy Elite ES-11S-C1 <sup>1</sup>	Group 1: Standard (K0816)
Weight Capacity	300 lbs.
Turning Radius	24.75"
Overall Width <sup>3</sup>	22.75"
Overall Length <sup>3</sup>	34.75"
Maximum Speed <sup>2</sup>	Up to 4 mph
Ground Clearance <sup>3</sup>	.75" at motors
Front Wheels	5" solid anti-tips
Drive Wheels	9" solid
Rear Wheels	6" solid casters
Drivetrain	Two-motor, in-line, front-wheel drive
Braking System	Regenerative and electro-mechanical
Suspension Type	Limited
High-back Seat Sizes	16" x 16"-18", 18" x 18"-20", 20" x 18"- 20"
Seat-to-Floor Range	21.5" - 23.5"
Per Charge Range <sup>3,5</sup>	Up to 12 miles
Standard Electronics	40A, PG GC3
Component Weights	Base 69 lbs. (heaviest piece) Seat 50.2 lbs. (comfort high-back) Batteries 25 lbs. each
Battery Requirements <sup>4,6</sup>	(2) 12 volt, U-1 batteries
Battery Charger	3.5 amp, off-board
Battery HCPCS	E2365
Warranty	Frame: 5-year limited Drivetrain: 13-month Electronics: 13-month Battery: 13-month

- 1 Jazzy® Power Chairs are FDA Class II Medical Devices designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 4 AGM or gel cell type required.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Battery weight may vary based on manufacturer.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

**PRIDE**  
MOBILITY  
LIVE YOUR BEST®

(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

0003661-102021

**Jazzy**<sup>® 1</sup>  
ELITE ES-1

**COMFORT AND  
PERFORMANCE  
FOR YOUR ACTIVE  
LIFESTYLE**



FDA Class II  
Medical Device<sup>1</sup>

**AMERICA'S FAVORITE POWER CHAIR**

gojazzy.com |    

# EFFORTLESS CONTROL AND COMFORT



The sleek and elegant **Jazzy® Elite ES-1<sup>1</sup>** delivers high-class maneuvering with front-wheel drive technology. Enjoy effortless controls and a comfortable, high-back seat on this power chair. Choose Jazzy every time.

- Depth-adjustable, limited-recline, fold-flat, high-back seat
- Universal, fixed, in-line joystick mount
- Flip-up foot platform
- Black, flat-free, non-marking tires
- Easy-access free-wheel levers
- 60" lap belt

## Colors



## Seating

Black, vinyl, depth-adjustable and limited recline standard comfort high-back with 60" lap belt available sizes:

- 16"W x 16" - 18"D
- 18"W x 18" - 20"D
- 20"W x 18" - 20"D

## Optional Seating

- Black vinyl, high-back seat up to 20" x 20"

## Options

- Elevating leg rests
- Heel loops
- Residual limb support
- Swing-away legrests
- Swing-away joystick

## Accessories

- Cane/crutch holder
- Cell phone holder
- Cup holder
- Flag holder
- Lap belt (50", 70", 80", 90")
- Oxygen tank holder
- Rear basket
- Saddle bag
- Walker holder
- Weather cover
- XLR USB charger



Front-wheel drive design

Easy-access freewheel levers

PG GC3 Controller

## Military Patches

