

VivaLIFT!®

POWER RECLINER COLLECTION¹

THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

*Wireless Charging Station**



Shown in Badlands Mushroom Fabric

WITH
COUNTLESS
POSITION
COMBINATIONS

ATLAS PLUS¹

When you need a moment for yourself, the **Atlas Plus Collection** is there for you. Customize your comfort while you watch the big game or read your favorite book. Enjoy a refreshing beverage and keep it close with the cup holder, and conveniently charge your phone with the wireless charger. Take advantage of those rare moments of solitude with the Atlas Plus Recliner.

* The wireless charger is only compatible with smart phones that have wireless charging capabilities. Additionally, no other smart devices, (i.e. smart watches) should be used with the wireless charger.



Scan here for
more information
on **VivaLift!**®



PRIDE®
MOBILITY
LIVE YOUR BEST®

ATLAS PLUS¹

VivaLIFT![®]
POWER RECLINER COLLECTION¹

COLLECTION	ATLAS PLUS ¹
MODEL	PLR-2985M
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift 5'4"-5'9"
SIZE	Medium
WEIGHT CAPACITY	400 lbs.
SEAT-TO-FLOOR	19.5"
SEAT DEPTH	21"
SEAT WIDTH	22"
TOP OF BACK TO SEAT	26"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	16"
OVERALL WIDTH	39"
CHAIR WEIGHT	164 lbs.
WARRANTY	7-Yr. Pro-rated

INDEPENDENT FUNCTIONS

-  1. The power headrest lifts your head, neck and shoulders for maximum support.
-  2. The upward lift offers convenient stand assist for getting in and out of the recliner.
-  3. The power lumbar enables you to personalize your comfort.
-  4. The footrest extension allows you to relax and stretch out completely.
-  5. The backrest moves up and down independently to help you find the perfect position.
-  6. Infinite position recline enables the backrest and footrest to move independently from one another so that a variety of positions can be reached.



STANDARD FEATURES

7. The USB remote adjusts all functions independently and allows you to charge your smart device.
8. Lithium battery backup ensures your recliner is operational when the power goes out.

1. VivaLift![®] Power Recliners are FDA Class II Medical Devices designed to aid individuals with mobility impairments.

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.

 = Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.

STANDARD FABRIC: BADLANDS



MUSHROOM

WALNUT

STEEL

(US) 800-800-8586 / (CAN) 888-570-1113
prideliftchairs.com / pridemobility.com

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]