

**Quantum® Rehab**

401 York Ave., Duryea, PA 18642

Phone: 866-800-2002 | Fax: 866-707-3422 | Email: [quantumorders@pridemobility.com](mailto:quantumorders@pridemobility.com)**Q6 Edge® 3 Order Form with Contoured, Synergy®, and TRU-Balance® Flex Seating**

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_

Provider Name: \_\_\_\_\_

Contact: \_\_\_\_\_

Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

Email: \_\_\_\_\_

PO Number: \_\_\_\_\_

Marked for: \_\_\_\_\_

Ship to Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_



TRU-Balance® 3 Power Lift Seating Shown

**INTRODUCTION**

This form is interactive when viewed with [Adobe Acrobat Reader](#) and may not function correctly if opened with applications other than Acrobat. Complete the form by placing checks in the desired boxes and provide information in the interactive fields. Buttons shown at the bottom of the form may be utilized to print or submit the order form through a desktop email application. To email via a web-based application, please 'Save As' and attach the PDF to your email.

This order form contains a large variety of options to fulfill various patient needs. Descriptions and section notes, such as optional or required, are included to help you complete your order. Please contact Quantum Sales at 866-800-2002 if assistance is needed. Send the completed order form by fax (866-707-3422) or email ([quantumorders@pridemobility.com](mailto:quantumorders@pridemobility.com)). Incomplete forms may delay the quote or order. Customer service will contact you if the order is incomplete or if there are compatibility issues. If special order requests are needed, be sure that the Patient Information section is completed or include a completed [physical assessment form](#) with this order form.

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. Payer coding, coverage, and bundling guidelines may apply. All prices are MSRP. Prices, specifications, part numbers, and availability are subject to change without notice. Prices and part numbers as shown are only available when configured as a complete power chair. Please contact technical service for accurate parts ordering. Options noted with "XRef" have multiple possible part numbers based on system configuration. Please see [cross-reference spreadsheet](#) for the XRef part numbers.


**Thank you for choosing Quantum!****PATIENT INFORMATION**

Weight: \_\_\_\_\_ Height: \_\_\_\_\_

A. Shoulder Width: \_\_\_\_\_ D. Max Sitting Height: \_\_\_\_\_ G. Thigh Depth: \_\_\_\_\_

B. Chest Width: \_\_\_\_\_ E. Shoulder Height: \_\_\_\_\_ H. Lower Leg Length: \_\_\_\_\_

C. Hip Width: \_\_\_\_\_ F. Axilla Height: \_\_\_\_\_ I. Elbow Height: \_\_\_\_\_

**1. BASE MODEL**

**REQUIRED.** Select a group 3 power base. 300 lb. weight capacity with standard 6 mph programmable motors and 8 amp off board charger. XLR USB charger is included. Dimensions of base without legrests: 24"Wx35.5"L

☐ Quantum Q6 Edge 3 3S-SS ..... \$7,895

Part: Q6 EDGE 3 3S-SS. HCPCS: K0848

Solid seat pan standard with no power option.

☐ Quantum Q6 Edge 3 3SP-SS ..... \$7,895

Part: Q6 EDGE 3 3SP-SS. HCPCS: K0856

Single power option for power elevating seat.

☐ Quantum Q6 Edge 3 3S-C ..... \$7,895

Part: Q6 EDGE 3 3S-C. HCPCS: K0849

Contoured seat (captain) standard.

☐ Quantum Q6 Edge 3 3SP-C ..... \$7,895

Part: Q6 EDGE 3 3SP-C. HCPCS: K0857

Contoured seat (captain), single power option for 12" lift.

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Marked for: \_\_\_\_\_

## 2. PATIENT'S WEIGHT

**REQUIRED.** The patient's weight is required to quote and order the power chair.

Patient's Weight: \_\_\_\_\_

## 3. BASE COLOR

**REQUIRED.** Select one shroud color. Colors shown here may differ from actual product depending on monitor/printer calibration used.

### 3.1. Match Colors

Match All Color Options to Selection Below

Part: N/A. Any subsequent color selections, if available, will be matched to color chosen below when the order is placed. Please note that some color selections, such as joystick armor, may have a MSRP.

### 3.2a. Standard Shrouds

|                                                                                                                                 |                                                                                                                                    |                                                                                                                                  |                                                                                                                                    |
|---------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------|
|  <b>Red Corvette</b><br>Part: DGN1801101-PRD1  |  <b>White Satin</b><br>Part: DGN1801101-PWH2      |  <b>Midnight Blue</b><br>Part: DGN1801101-PBU10 |  <b>Military Green</b><br>Part: DGN1801101-PGN8 |
|  <b>Cherry Bomb</b><br>Part: DGN1801101-PRD4   |  <b>Deep Purple</b><br>Part: DGN1801101-PPR5      |  <b>Electric Blue</b><br>Part: DGN1801101-PBU9  |  <b>Sugar Plum</b><br>Part: DGN1801101-PPR7     |
|  <b>Back in Black</b><br>Part: DGN1801101-PBK1 |  <b>Raspberry Beret</b><br>Part: DGN1801101-PPK4  |  <b>Blue Jean</b><br>Part: DGN1801101-PBU11     | <b>To Be Determined</b><br>Follow-up required.                                                                                     |
|  <b>Night Train</b><br>Part: DGN1801101-PBK3   |  <b>Pretty in Pink</b><br>Part: DGN1801101-PPK3   |  <b>Lemon Crush</b><br>Part: DGN1801101-PYL4    |                                                                                                                                    |
|  <b>Grey Street</b><br>Part: DGN1801101-PSL3   |  <b>Evolution Orange</b><br>Part: DGN1801101-POR7 |  <b>Green Machine</b><br>Part: DGN1801101-PGN11 |                                                                                                                                    |

### 3.2b. Hydrographic Shrouds

|                                                                                                                                                                              |                                                                                                                                                                    |                                                                                                                                                                       |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/>  <b>Military Digital Camo</b> ..... \$250<br>Part: DGN1801101-WCM7 | <input type="checkbox"/>  <b>Faded Glory</b> ..... \$250<br>Part: DGN1801101-WFL2 | <input type="checkbox"/>  <b>Carbon Fiber</b> ..... \$250<br>Part: DGN1801101-WCB1 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## 4. FENDER LIGHTS

**Optional.** Fender lights operate by pressing included single switch button.

☐ **Switch Fender Lights** ..... No Charge  
Part: DGN154117

## 5. FULL LIGHT PACKAGES

**Optional.** Requires NE+ 8 Key or Q-Logic 3 Joystick for Full Light Package. Select one light package below.

- ☐ **Base Mounted, LED, Full Light Package** ..... \$995  
Part without iLevel: ELE1705999. Part with iLevel: ELE168643  
Fender lights will be omitted with Base Mounted, LED, Full Light Package.
- ☐ **Seat Mounted, LED, Full Light Package** ..... \$995  
Flex part: ELE131477. Synergy part: ELEASMB5719  
Not available with contoured seating.



## 6. DRIVE MOTORS

**REQUIRED.**

|                                                             |                                                                                                                                                                                                       |
|-------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>6 mph Drive Motors</b> ..... Standard<br>Part: MOT147433 | <b>Accu-Trac Motor Technology</b> ..... \$1,295<br>Part: MOT147881. HCPCS: K0108<br>Requires Q-Logic 3 electronics. Delivers enhanced traction and tracking, especially with switched drive controls. |
|-------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

## 7. BASE MOUNTED FOOT PLATFORM

**REQUIRED.** May interfere with seat mounted legrests. Not available with center mount foot platforms.

|                                                                                    |                                                                             |
|------------------------------------------------------------------------------------|-----------------------------------------------------------------------------|
| <b>Base Mounted Foot Platform</b> ..... Standard<br>Part: FRMASMB16077             | <b>Omit base mounted foot platform</b> ..... No Charge<br>Part: PTOINDV3674 |
| <b>Bell Shaped, Base Mounted Foot Platform</b> ..... \$235<br>Part: RIGASMB7120041 |                                                                             |

## 8. BATTERIES

**REQUIRED.** The powerchair requires 2 batteries and will be installed if ordered.

- ☐ **NF-22, Gel Cell, Interceptor Battery** ..... \$310 ea.  
Part: BAT1704613. HCPCS: E2361
- ☐ **NF-22, AGM, MK Battery** ..... \$395 ea.  
Part: BAT121346. HCPCS: E2361
- ☐ **No batteries.** Provider will assume responsibility of batteries.  
Part: Configuration selection
- ☐ **NF-22, Gel Cell, MK Battery** ..... \$395 ea.  
Part: BATGEL1008. HCPCS: E2361

**9. DRIVE WHEELS****REQUIRED.** Select a type of drive wheels and accent color. Knobby tires increase overall base width by 1.5".**9.1. Drive Wheels**

|                                                                               |                                                                               |
|-------------------------------------------------------------------------------|-------------------------------------------------------------------------------|
| 3"x14" Flat Free Black Tires ..... Standard<br>Part: PTOINDV3525 (pair)       | 3"x14" Pneumatic Black Tires ..... No Charge<br>Part: WHL1713925 x2           |
| 4"x14" Knobby, Flat Free Black Tires ..... \$275.50 ea.<br>Part: WHL147601 x2 | 4"x14" Knobby, Pneumatic Black Tires ..... \$275.50 ea.<br>Part: WHL147599 x2 |

**9.2. Drive Wheel Accent Color**

Not available with knobby tires.

|                                                         |                                                             |                                                                 |
|---------------------------------------------------------|-------------------------------------------------------------|-----------------------------------------------------------------|
| Back in Black ..... Standard<br>Part: DGN164245-PBK1 x2 | Raspberry Beret ..... No Charge<br>Part: DGN164245-PPK4 x2  | Green Machine ..... No Charge<br>Part: DGN164245-PGN11 x2       |
| Red Corvette ..... No Charge<br>Part: DGN164245-PRD1 x2 | Pretty in Pink ..... No Charge<br>Part: DGN164245-PPK3 x2   | Military Green ..... No Charge<br>Part: DGN164245-PGN8 x2       |
| Cherry Bomb ..... No Charge<br>Part: DGN164245-PRD4 x2  | Evolution Orange ..... No Charge<br>Part: DGN164245-POR7 x2 | Sugar Plum ..... No Charge<br>Part: DGN164245-PPR7 x2           |
| Night Train ..... No Charge<br>Part: DGN164245-PBK3 x2  | Midnight Blue ..... No Charge<br>Part: DGN164245-PBU10 x2   | Military Digital Camo ..... \$50 ea.<br>Part: DGN164245-WCM7 x2 |
| Grey Street ..... No Charge<br>Part: DGN164245-PSL3 x2  | Electric Blue ..... No Charge<br>Part: DGN164245-PBU9 x2    | Faded Glory ..... \$50 ea.<br>Part: DGN164245-WFL2 x2           |
| White Satin ..... No Charge<br>Part: DGN164245-PWH2 x2  | Blue Jean ..... No Charge<br>Part: DGN164245-PBU11 x2       | Carbon Fiber ..... \$50 ea.<br>Part: DGN164245-WCB1 x2          |
| Deep Purple ..... No Charge<br>Part: DGN164245-PPR5 x2  | Lemon Crush ..... No Charge<br>Part: DGN164245-PYL4 x2      | To Be Determined<br>Follow-up required.                         |

**10. SEAT TOWERS/POWER ADJUSTABLE SEAT HEIGHT****REQUIRED.** Select one option below.**10.1a. Power Base Seat Towers**

|                                                                                      |                                                                                        |
|--------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| Standard Seat Towers ..... Standard<br>Part: Configuration selection, no part number | Extended Seat Towers (increased seat to floor height) ... No Charge<br>Part: MEC120301 |
|--------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|

**10.1b. Power Adjustable Seat Height**

The seat to floor height will be increased. Lift/elevate will be operated through the joystick. Power elevating seat includes a 1/4 speed inhibit.

|                                                                                                                                                                                                                                                                                                                           |                                                                                                                                                                                                               |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 12" Power Lift for Captain Seat with iLevel ..... \$2,600<br>Part: PTOINDV3857. HCPCS: E2300<br>For use with K0857 base code. Must select <u>Power Lift Captain Seat</u> in contoured seat section. Includes armrests and lap belt. Available with swing away legrests or manually adjustable center mount foot platform. | Power Elevating Seat with 7" Stroke ..... \$1,895<br>Contoured seat part: FRMASMB16310. HCPCS: E2300<br>Comfort contoured seat part: MEC133630. HCPCS: E2300<br>Synergy seat part: FRMASMB16312. HCPCS: E2300 |
|                                                                                                                                                                                                                                                                                                                           | Power Elevating Seat with 6" Stroke ..... \$1,895<br>Contoured seat part: FRMASMB16309. HCPCS: E2300<br>Comfort contoured seat part: MEC133629. HCPCS: E2300<br>Synergy seat part: FRMASMB16311. HCPCS: E2300 |

**11. ELECTRONICS****REQUIRED.** Select a type of drive control from the following sections. The No Joystick option is used to omit the joystick and utilize another control, such as a specialty control at the end of this form or found on the [Alternative Controls and Electronics Order Form](#).**11.1a. Non-Expandable Controls**

Select a joystick below or omit selection if another controller type is desired. Non-expandable does not allow supplementary drive controllers.

|                                                                                                                                                                                   |                                                                                                                                                                |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> NE 4-Key Joystick ..... Standard<br>Part: CTLDC1554<br>Not available with power seating.                                                                 | <input type="checkbox"/> NE+ 6-Key Joystick (2 Actuators) ..... \$200<br>Part: CTLDC1560<br>Control up to 2 actuators through joystick.                        |
| <input type="checkbox"/> NE+ 6-Key Joystick (1 Actuator) ..... No Charge<br>Part: CTLDC1560<br>Control of one actuator through joystick on single power base model (3SP-SS) only. | <input type="checkbox"/> NE+ 8-Key Joystick (2 Actuators) ..... \$400<br>Part: CTLDC1561<br>Control up to 2 actuators through joystick and full light package. |
|                                                                                                                                                                                   | <input type="checkbox"/> No Joystick, Setup for NE+ Joystick ..... No Charge<br>Part: Configuration selection                                                  |

**11.1b. Expandable Controls - Quantum Intelligent Systems™**

Select expandable controls or omit selection if another controller type is desired. Expandable controls allow multiple or specialized controllers.

**1. Q-Logic 3 Required Components**

Q-Logic 3 expandable controls require the controller module and expandable harness below. If a Q-Logic joystick or specialty control is selected, Customer Service will quote the order with these items.

|                                                                                                     |                                                                                                               |
|-----------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> Q-Logic 3 EX Controller Module ..... \$695<br>Part: See XRef. HCPCS: E2377 | <input type="checkbox"/> Q-Logic 3 Expandable Controller Harness ..... \$425<br>Part: HAR144005. HCPCS: E2313 |
|-----------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|

**2. Q-Logic 3 Joystick**

|                                                                                           |                                                                                               |
|-------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------|
| <input type="checkbox"/> Q-Logic 3 EX Joystick ..... \$700<br>Part: CTL167724             | <input type="checkbox"/> Q-Logic 3 EX Joystick for Full Lights ..... \$700<br>Part: CTL167725 |
| <input type="checkbox"/> No Joystick, Setup for Q-Logic<br>Part: Configuration selection. |                                                                                               |



Electronics section continued on following page.

**11. ELECTRONICS CONTINUED****11.2. Q-Logic 3 Joystick Armor****REQUIRED with Q-Logic 3 Joystick selections above.** Reference colors shown on page 2. \*Night train joystick armor will be the black, soft touch.

**Black, Soft Touch** ..... Standard  
 Part: DGN124307  
 Also matches Night Train shrouds.

**Non-matching Armor** ..... \$50  
 Part: PTOINDV3422  
 Select color below

**Match Armor to Base Color** ..... \$50  
 Part: PTOINDV3421  
 Part below will be selected with system configurator.

|                                                                        |                                                |                                                 |                                               |                                                      |
|------------------------------------------------------------------------|------------------------------------------------|-------------------------------------------------|-----------------------------------------------|------------------------------------------------------|
| <b>Red Corvette</b><br>Part: DGN124336-PRD1                            | <b>White Satin</b><br>Part: DGN124336-PWH2     | <b>Evolution Orange</b><br>Part: DGN124336-POR7 | <b>Lemon Crush</b><br>Part: DGN124336-PYL4    | <b>Military Digital Camo</b><br>Part: DGN124336-WCM7 |
| <b>Cherry Bomb</b><br>Part: DGN124336-PRD4                             | <b>Deep Purple</b><br>Part: DGN124336-PPR5     | <b>Midnight Blue</b><br>Part: DGN124336-PBU10   | <b>Green Machine</b><br>Part: DGN124336-PGN11 | <b>Faded Glory</b><br>Part: DGN124336-WFL2           |
| <b>Back in Black</b><br>Part: DGN124336-PBK1                           | <b>Raspberry Beret</b><br>Part: DGN124336-PPK4 | <b>Electric Blue</b><br>Part: DGN124336-PBU9    | <b>Military Green</b><br>Part: DGN124336-PGN8 | <b>Carbon Fiber</b><br>Part: DGN124336-WCB1          |
| <input type="checkbox"/> <b>Night Train*</b><br>*See Black, Soft Touch | <b>Pretty in Pink</b><br>Part: DGN124336-PPK3  | <b>Blue Jean</b><br>Part: DGN124336-PBU11       | <b>Sugar Plum</b><br>Part: DGN124336-PPR7     | <b>To Be Determined</b><br>Follow-up required.       |
| <b>Grey Street</b><br>Part: DGN124336-PSL3                             |                                                |                                                 |                                               |                                                      |

**12. JOYSTICK HANDLES****Optional.** The NE/NE Plus Joysticks are standard with a carrot handle (See A in image). The optional handle will be mounted to the joystick and the standard handle will be shipped with the unit.


- |                                                                                                                      |                                                                                                             |
|----------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <b>Mushroom Handle (B)</b> ..... \$85<br>Part: ACCASMB1860. HCPCS: E2323                    | <input type="checkbox"/> <b>Extended I-Handle (G)</b> ..... \$85<br>Part: ACCASMB1857. HCPCS: E2323         |
| <input type="checkbox"/> <b>Chin Cup (C)</b> ..... \$85<br>Part: ACCASMB1861. HCPCS: E2323 (E2324 when chin control) | <input type="checkbox"/> <b>Bodypoint 3" Goal Post Handle</b> ..... \$99<br>Part: PLSKNCP1057. HCPCS: E2323 |
| <input type="checkbox"/> <b>T-Bar (D)</b> ..... \$120<br>Part: ACCASMB1858. HCPCS: E2323                             | <input type="checkbox"/> <b>Bodypoint 4" Goal Post Handle</b> ..... \$99<br>Part: PLSKNCP1058. HCPCS: E2323 |
| <input type="checkbox"/> <b>Medium Ball Handle (E)</b> ..... \$65<br>Part: ACCASMB1859. HCPCS: E2323                 | <input type="checkbox"/> <b>Q-Logic 3 Carrot Handle (A)</b> ..... No Charge<br>Part: DGN124308              |
| <input type="checkbox"/> <b>Soft Ball Handle (F)</b> ..... \$65<br>Part: ACCASMB1856. HCPCS: E2323                   |                                                                                                             |

**13. JOYSTICK MOUNTING BRACKETS****Optional.** Select a mount to be used with a joystick selected in previous section.

Omit Joystick Mounting Bracket


**13.1a. Fixed, In-line Joystick Bracket**

- ☐ **Fixed, In-line, Bracket** ..... Standard  
 Q-Logic Part: FRMASMB10617. NE/NE+ Part: FRMASMB12744. SAJ Part: FRMASMB10454
- ☐ **Left Side Mount** .....  
 Part: PTOINDV1225
- ☐ **Right Side Mount** .....  
 Part: PTOINDV1226

**13.1b. Swing-Away (Retractable) Joystick Bracket**


Allows the joystick to retract to the side of the armrest while maintaining the same joystick orientation. "Inline bracket" will be offset 1" when used with tube style arms. Select "Offset Bracket" to have the joystick inline with armpad on tube style armrests.

- ☐ **Swing-away, Inline, Left Side Mount** ..... \$245  
 Part: See XRef. HCPCS: E1028
- ☐ **Swing-away, 1" Offset, Left Side Mount** ..... \$245  
 Part: See XRef. HCPCS: E1028
- ☐ **Swing-away, Inline, Right Side Mount** ..... \$245  
 Part: See XRef. HCPCS: E1028
- ☐ **Swing-away, 1" Offset, Right Side Mount** ..... \$245  
 Part: See XRef. HCPCS: E1028

**13.1c. Specialty Joystick Mounts**

Select a bracket and side to mount it.

- ☐ **Retract4 Mount** ..... \$350  
 NE/NE+ Part: MEC138876. HCPCS: E1028  
 Q-Logic Part: MEC138860. HCPCS: E1028  
 Full 180° of movement that translates into multi-position use, straight, right offset, left offset, or fold away. Joystick orientation varies as bracket is adjusted.
- ☐ **Stealth Swing-Away/Flip Down Bracket** ..... \$294.25  
 NE/NE+ Part: ST-TWBM480QL-Q. HCPCS: E1028  
 Q-Logic Part: ST-TWBM480QL-Q. HCPCS: E1028  
 SAJ Part: ST-TWBM480CJ-Q. HCPCS: E1028  
 Allows the joystick to pivot to side of armrest, or flip-down if bracket is adjusted.
- Multi-axis Mount** ..... \$100  
 NE/NE+ Part: ACC142528. HCPCS: E1028  
 Q-Logic Part: ACC131216. HCPCS: E1028  
 SAJ Part: ACC142740. HCPCS: E1028  
 Allows height adjustment of joystick.

**Specialty Mount Side**

- ☐ **Left Side Mount** .....  
 Part: PTOINDV1225
- ☐ **Right Side Mount** .....  
 Part: PTOINDV1226

Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Marked for: \_\_\_\_\_

#### 14. CONTOURED SEATING

**Optional.** Select only one seating system. Special order seat sizes can be requested by checking the seat type and writing in the seat size in the special order seat size.

##### 14.1. Occupied Transit, WC-19

WC-19 Occupied Transit Option for Contoured Seat ..... \$350  
Part: KIT122609

Only available with Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline. To determine if your chair is WC-19 compliant (occupied), please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>

##### 14.2a. Captain Style Seating

###### 1a. Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline

Comfort seat, full length armrests must be selected with a depth adjustable high back.

|                                                              |                                                              |                                                              |
|--------------------------------------------------------------|--------------------------------------------------------------|--------------------------------------------------------------|
| 16"Wx16"-18"D w/ 60" Lap Belt ..... \$150<br>Part: SET131843 | 18"Wx18"-20"D w/ 60" Lap Belt ..... \$150<br>Part: SET131844 | 20"Wx18"-20"D w/ 60" Lap Belt ..... \$150<br>Part: SET131845 |
|--------------------------------------------------------------|--------------------------------------------------------------|--------------------------------------------------------------|

###### 1b. Depth Adjustable High Back with Limited Recline

Full length armrests must be selected with a depth adjustable high back.

|                                                                      |                                                                  |
|----------------------------------------------------------------------|------------------------------------------------------------------|
| 18"Wx18"-20"D with 60" Lap Belt ..... No Charge<br>Part: SETASMB6895 | 22"Wx20"-22"D with 70" Lap Belt ..... \$150<br>Part: SETASMB6897 |
| 20"Wx18"-20"D with 60" Lap Belt ..... \$150<br>Part: SETASMB6896     |                                                                  |

###### 1c. Power Lift, Fold Flat Captain Seat

Only available with 12" Power Lift for Captain Seat with iLevel. Includes armrests and 60" lap belt.

|                                               |                                               |                                               |
|-----------------------------------------------|-----------------------------------------------|-----------------------------------------------|
| 16"Wx16"-18"D ..... \$150<br>Part: SET1701577 | 18"Wx18"-20"D ..... \$150<br>Part: SET1701578 | 20"Wx20"-22"D ..... \$150<br>Part: SET1701579 |
|-----------------------------------------------|-----------------------------------------------|-----------------------------------------------|

###### 1d. High Back with Limited Recline

|                                            |                                            |                                            |
|--------------------------------------------|--------------------------------------------|--------------------------------------------|
| 16"Wx16"D ..... \$150<br>Part: SETASMB2770 | 16"Wx22"D ..... \$150<br>Part: SETASMB4887 | 20"Wx16"D ..... \$150<br>Part: SETASMB2780 |
| 16"Wx18"D ..... \$150<br>Part: SETASMB2771 | 18"Wx16"D ..... \$150<br>Part: SETASMB2775 | 20"Wx22"D ..... \$150<br>Part: SETASMB4889 |
| 16"Wx20"D ..... \$150<br>Part: SETASMB2772 | 18"Wx22"D ..... \$150<br>Part: SETASMB4888 |                                            |

##### 14.2b. Contoured Seat with Solid Seat Pan

###### 1a. Comfort, Depth Adjustable, Height Adjustable, High Back with Limited Recline and Solid Seat Pan

Comfort seat, full length armrests must be selected with a depth adjustable high back.

|                                                                |                                                                |
|----------------------------------------------------------------|----------------------------------------------------------------|
| 16"Wx16"-18"D with 60" Lap Belt ..... \$150<br>Part: SET145814 | 18"Wx18"-20"D with 60" Lap Belt ..... \$150<br>Part: SET145815 |
|----------------------------------------------------------------|----------------------------------------------------------------|

###### 1b. Power Lift, Fold Flat Contoured Seat with Solid Seat Pan

Only available with 12" Power Lift with iLevel. Includes armrests and lap belt.

|                                                               |                                                               |                                                               |
|---------------------------------------------------------------|---------------------------------------------------------------|---------------------------------------------------------------|
| 16"Wx16"-18"D w/ 60" Lap Belt ..... \$150<br>Part: SET1701580 | 18"Wx18"-20"D w/ 60" Lap Belt ..... \$150<br>Part: SET1701581 | 20"Wx20"-22"D w/ 60" Lap Belt ..... \$150<br>Part: SET1701582 |
|---------------------------------------------------------------|---------------------------------------------------------------|---------------------------------------------------------------|

###### 1c. High Back with Limited Recline and Solid Seat Pan

|                                            |                                            |                                            |
|--------------------------------------------|--------------------------------------------|--------------------------------------------|
| 16"Wx16"D ..... \$150<br>Part: SETASMB3178 | 18"Wx16"D ..... \$150<br>Part: SETASMB3181 | 20"Wx16"D ..... \$150<br>Part: SETASMB3183 |
| 16"Wx18"D ..... \$150<br>Part: SETASMB3179 | 18"Wx18"D ..... \$150<br>Part: SETASMB2848 | 20"Wx18"D ..... \$150<br>Part: SETASMB3184 |
| 16"Wx20"D ..... \$150<br>Part: SETASMB3180 | 18"Wx20"D ..... \$150<br>Part: SETASMB3182 | 20"Wx20"D ..... \$150<br>Part: SETASMB3185 |
| 16"Wx22"D ..... \$150<br>Part: SETASMB5006 | 18"Wx22"D ..... \$150<br>Part: SETASMB5007 | 20"Wx22"D ..... \$150<br>Part: SETASMB5008 |

##### 14.3. Seat to Floor Height

Required with any contoured seating selection.

|                                           |                                                                                         |                                            |
|-------------------------------------------|-----------------------------------------------------------------------------------------|--------------------------------------------|
| Lowest Seat to Floor<br>Part: PTOINDV3558 | Middle Seat to Floor<br>Part: PTOINDV3559<br>Not available with Comfort high back seat. | Highest Seat to Floor<br>Part: PTOINDV3560 |
|-------------------------------------------|-----------------------------------------------------------------------------------------|--------------------------------------------|

##### 14.4. Set Depth Selection for Comfort Seat

Required with Comfort Seat selection above.

|                                |                                |                                |                                |
|--------------------------------|--------------------------------|--------------------------------|--------------------------------|
| 14" Depth<br>Part: SETDPTH1004 | 16" Depth<br>Part: SETDPTH1006 | 18" Depth<br>Part: SETDPTH1008 | 20" Depth<br>Part: SETDPTH1010 |
|--------------------------------|--------------------------------|--------------------------------|--------------------------------|

##### 14.5. Back Height Selection for Comfort Seat

Required with Comfort Seat selection above.

|                                                         |                                                          |
|---------------------------------------------------------|----------------------------------------------------------|
| Set Back Height to Lowest Position<br>Part: SETHGHT1100 | Set Back Height to Highest Position<br>Part: SETHGHT1101 |
|---------------------------------------------------------|----------------------------------------------------------|

##### 14.6. Special Order Contoured Seat

Complete information below. Contact inside sales specialist for pricing.

Type of Seat: \_\_\_\_\_ Width: \_\_\_\_\_ Depth: \_\_\_\_\_

Contoured seating section continued on following page.

**14. CONTOURED SEATING CONTINUED****14.7. Contoured Seat Flip-Up, Height Adjustable, Removable Armrests****REQUIRED except with fold flat, power lift captain seat.** Select a left and right armrest or to omit them.

|                                                                         |                                                                                        |
|-------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| Full Length, Left Armrest.....\$85<br>Part: FRMASMB8614. HCPCS: E0973   | Comfort Seat, Right, Full Length Armrest.....\$85<br>Part: FRMASMB16315. HCPCS: E0973  |
| Desk Length, Left Armrest.....\$85<br>Part: FRMASMB8612. HCPCS: E0973   | Comfort Seat, Left, Full Length, Armrest .....\$85<br>Part: FRMASMB16316. HCPCS: E0973 |
| Full Length, Right Armrest.....\$85<br>Part: FRMASMB8615. HCPCS: E0973  | No Armrest ..... Standard<br>Part: N/A                                                 |
| Desk Length, Right Armrest .....\$85<br>Part: FRMASMB8613. HCPCS: E0973 |                                                                                        |

**15. SYNERGY SEATING SYSTEM****Optional.** Widths, depths, and heights are available in 1" increments in the noted ranges. Maximum seat size 22"W x 22"D.**15.1. Synergy Frame Selection**

|                              |       |
|------------------------------|-------|
| Synergy Seating System ..... | \$434 |
| Part: PTOINDV2995            |       |

**15.2. Synergy Transit Kit****Optional.** The power base is standard with Unoccupied Transit Loops. 60" Sure-Lok Lap Belt is required with Occupied Transit Option.

|                                                                              |       |
|------------------------------------------------------------------------------|-------|
| <input type="checkbox"/> WC-19 Occupied Transit Option for Synergy Seat..... | \$275 |
| Part: KITASMB1830                                                            |       |

To determine if your chair is WC-19 compliant, please refer to this website: <http://wc-transportation-safety.umtri.umich.edu/crash-tested-product-lists/wheelchairs>**15.3. Seat Dimension - Width****REQUIRED with Synergy seating.** Part number is SETWDTH plus number following # below.

|                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| 10" W<br>(#1000) | 11" W<br>(#1001) | 12" W<br>(#1002) | 13" W<br>(#1003) | 14" W<br>(#1004) | 15" W<br>(#1005) | 16" W<br>(#1006) | 17" W<br>(#1007) | 18" W<br>(#1008) | 19" W<br>(#1009) | 20" W<br>(#1010) | 21" W<br>(#1011) | 22" W<br>(#1012) |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|

**15.4. Seat Dimension - Width Package****REQUIRED with Synergy seating.** Maximum seat width is 22".

|                                             |                                             |                                             |                                                                    |
|---------------------------------------------|---------------------------------------------|---------------------------------------------|--------------------------------------------------------------------|
| 10"-14" Width Package<br>Part: FRMASMB18011 | 14"-18" Width Package<br>Part: FRMASMB18012 | 16"-20" Width Package<br>Part: FRMASMB18013 | 20"-24" Width Package<br>Part: FRMASMB18014 (Maximum width is 22") |
|---------------------------------------------|---------------------------------------------|---------------------------------------------|--------------------------------------------------------------------|

**15.5. Seat Dimension - Depth****REQUIRED with Synergy seating.** Part number is SETDPTH plus number following # below.

|                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |                  |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| 10" D<br>(#1000) | 11" D<br>(#1001) | 12" D<br>(#1002) | 13" D<br>(#1003) | 14" D<br>(#1004) | 15" D<br>(#1005) | 16" D<br>(#1006) | 17" D<br>(#1007) | 18" D<br>(#1008) | 19" D<br>(#1009) | 20" D<br>(#1010) | 21" D<br>(#1011) | 22" D<br>(#1012) |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|

**15.6. Seat Dimension - Depth Package****REQUIRED with Synergy seating.** Maximum seat depth is 22".

|                                            |                                             |                                            |                                                                   |
|--------------------------------------------|---------------------------------------------|--------------------------------------------|-------------------------------------------------------------------|
| 10"-14" Depth Package<br>Part: FRMASMB6981 | 14"-18" Depth Package<br>Part: FRMASMB17823 | 15"-20" Depth Package<br>Part: FRMASMB7732 | 20"-24" Depth Package<br>Part: FRMASMB7733 (Maximum depth is 22") |
|--------------------------------------------|---------------------------------------------|--------------------------------------------|-------------------------------------------------------------------|

**15.7. Synergy Seat to Floor****REQUIRED with Synergy seating.**

|                              |                              |                              |                              |                              |                              |
|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|
| 15.875"<br>Part: PTOINDV3652 | 16.875"<br>Part: PTOINDV3653 | 17.875"<br>Part: PTOINDV3654 | 18.875"<br>Part: PTOINDV3761 | 19.875"<br>Part: PTOINDV3762 | 20.875"<br>Part: PTOINDV3763 |
|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|------------------------------|

**15.8. Synergy Back Options****REQUIRED with Synergy seating.** See the Synergy and TRU-Comfort back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time.**1. Back Cane Selection**

|                                                                                         |           |                                                                                                                                                                                                                                                       |       |
|-----------------------------------------------------------------------------------------|-----------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|
| Synergy Back Cane .....                                                                 | No Charge | 2" Endomorph (1" per side) Back Cane.....                                                                                                                                                                                                             | \$195 |
| Part: FRMASMB15845                                                                      |           | Part: FRMASMB12115                                                                                                                                                                                                                                    |       |
| TRU-Comfort Back Canes .....                                                            | No Charge | 4" Endomorph (2" per side) Back Cane.....                                                                                                                                                                                                             | \$195 |
| Part: N/A, Configuration selection                                                      |           | Part: FRMASMB12139                                                                                                                                                                                                                                    |       |
| Flip Forward Back Cane .....                                                            | \$150     | Manual Recline Back .....                                                                                                                                                                                                                             | \$750 |
| Part: FRMASMB15846                                                                      |           | Part: FRMASMB12910. HCPCS: E1226                                                                                                                                                                                                                      |       |
| Not available with a vent tray, oxygen holder, or cane mounted and cantilever armrests. |           | Select an overall back height. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, rear accessories, power elevating seat, or articulating foot platform. Pivot point is 2.5" above the seat pan. |       |
| 2" Mesomorph (1" per side) Back Cane.....                                               | \$195     |                                                                                                                                                                                                                                                       |       |
| Part: FRMASMB12067                                                                      |           |                                                                                                                                                                                                                                                       |       |
| 4" Mesomorph (2" per side) Back Cane.....                                               | \$195     |                                                                                                                                                                                                                                                       |       |
| Part: FRMASMB12091                                                                      |           |                                                                                                                                                                                                                                                       |       |

Synergy seating section continued on following pages.

**15. SYNERGY SEATING SYSTEM CONTINUED****2. Synergy Back Cane Height**

Omit this selection if a TRU-Comfort back is desired. 16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests. 14" and 15" back canes are not available without push handles.

14" High Canes  
Part: FRMASMB15832

15" High Canes  
Part: FRMASMB15833

16" High Canes  
Part: FRMASMB15834

17" High Canes  
Part: FRMASMB15835

18" High Canes  
Part: FRMASMB15836

19" High Canes  
Part: FRMASMB15837

20" High Canes  
Part: FRMASMB15838

21" High Canes  
Part: FRMASMB15839

22" High Canes  
Part: FRMASMB15840

23" High Canes  
Part: FRMASMB15841

24" High Canes  
Part: FRMASMB15842

25" High Canes  
Part: FRMASMB15843

**3. Back Selection**

Synergy Back..... \$275  
Part: See XRef

TRU-Comfort, 20" Width or less..... \$675  
Part: See XRef. HCPCS: E2620

TRU-Comfort, 22" Width or greater..... \$775  
Part: See XRef. HCPCS: E2621

Sling Back..... \$98  
Part: See XRef

Available 18"-25" overall heights. Back cane height must match back height. Not available with flip up cantilever arms.

Omit Back..... No Charge  
Part: PTOINDV1287

**4. Synergy Overall Back Height**

Parts: See XRef, configuration selection with back selection.

12" High    15" High    18" High    21" High    24" High

13" High    16" High    19" High    22" High    25" High

14" High    17" High    20" High    23" High

|                |     | Width |     |     |      |      |      |      |      |     |      |
|----------------|-----|-------|-----|-----|------|------|------|------|------|-----|------|
|                |     | 12"   | 14" | 15" | 16"  | 17"  | 18"  | 19"  | 20"  | 21" | 22"  |
| Overall Height | 15" | S     | S   | -   | -    | -    | -    | -    | -    | -   | -    |
|                | 16" | S     | S   | -   | -    | -    | -    | -    | -    | -   | -    |
|                | 17" | -     | S   | S   | S    | S    | S    | S    | S    | S   | S    |
|                | 18" | -     | S   | S   | S    | S    | S    | S    | S    | S   | S    |
|                | 19" | -     | S   | S   | S    | -    | S    | S    | S    | S   | S    |
|                | 20" | -     | S   | S   | S TC | -    | S TC | S    | S TC | S   | S    |
|                | 21" | -     | S   | S   | S TC | S    | S TC | S    | S TC | S   | S    |
|                | 22" | -     | S   | S   | S TC | S TC | S TC | S TC | S TC | S   | S    |
|                | 23" | -     | S   | S   | S TC | S TC | S TC | S TC | S TC | S   | S    |
|                | 24" | -     | S   | S   | S TC | S TC | S TC | S TC | S TC | S   | S TC |
|                | 25" | -     | S   | S   | S TC | S TC | S TC | S TC | S TC | S   | S TC |

S Synergy Back    S TC Synergy or TRU-Comfort Back

**Synergy and TRU-Comfort Back Availability Chart**

**15.9. Synergy Armrests**

**REQUIRED with Synergy seating.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.

**1a. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

☐ ☐ Full Length, Flip Back Armrests..... Standard  
L R Left Part: FRMASMB7038. Right part: FRMASMB7039

☐ ☐ Desk Length, Flip Back Armrests ..... Standard  
L R Left part: FRMASMB7040. Right part: FRMASMB7041

**1b. Two-Post, Flip Back, Height Adjustable, Removable Armrests**

Not available with a depth less than 16". Required with flip forward back and laptrays.

☐ ☐ Full Length, Two-Post Armrests.....\$150 ea.  
L R Left Part: FRMASMB7062. Right part: FRMASMB7063. HCPCS: E0973

☐ ☐ Desk Length, Two-Post Armrests .....\$150 ea.  
L R Left part: FRMASMB7060. Right part: FRMASMB7061. HCPCS: E0973

**1c. Single-Post, Heavy Duty, 8"-12" Height Adjustable, Removable Armrests**

☐ ☐ Full Length, Heavy Duty Armrests, 8"-12" H ..... \$75 ea.  
L R Left Part: FRMASMB9982. Right part: FRMASMB9983. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 8"-12"H..... \$75 ea.  
L R Left part: FRMASMB9984. Right part: FRMASMB9985. HCPCS: E0973

**1d. Single-Post, Heavy Duty, 10"-14" Height Adjustable, Removable Armrests**

☐ ☐ Full Length, Heavy Duty Armrests, 10"-14"H..... \$75 ea.  
L R Left Part: FRMASMB9986. Right part: FRMASMB9987. HCPCS: E0973

☐ ☐ Desk Length, Heavy Duty Armrests, 10"-14"H ..... \$75 ea.  
L R Left part: FRMASMB9988. Right part: FRMASMB9989. HCPCS: E0973

**1e. Pediatric, Removable, Height Adjustable Armrests**

Only available with desk length armrests.

☐ ☐ Pediatric, Removable Armrests ..... \$75 ea.  
L R Left Part: FRMASMB10608. Right part: FRMASMB10609. HCPCS: E0973

**1f. Height Adjustable, Cantilever Armrests**

Height adjustment from 11"-17" depends on cane height. Not available with manual recline, flip forward back canes, or cane clamp pushbar.

☐ ☐ Full Length, Cantilever Armrests ..... \$97.50 ea.  
L R Left Part: FRMASMB16188. Right part: FRMASMB16189

☐ ☐ Desk Length, Cantilever Armrests..... \$97.50 ea.  
L R Left part: FRMASMB16190. Right part: FRMASMB16191

**1g. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

☐ ☐ Full Length, Single Post Armrests ..... \$75 ea.  
L R Left Part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973

☐ ☐ Desk Length, Single Post Armrests ..... \$75 ea.  
L R Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973

**15.10. Synergy Armpads**

**REQUIRED with armrest selection above.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.

Straight Armpads..... Standard  
Desk/Desk part: UPSASMB1022. Full/Full part: UPSASMB1021.  
Desk/Full part: UPSASMB1023

Waterfall Armpads..... No Charge  
Desk/Desk part: UPSASMB1025. Full/Full part: UPSASMB1024.  
Desk/Full part: UPSASMB1026

No Armpad Selection ..... No Charge  
Part: PTOINDV1373

*Synergy seating section continued on following page.*

**15. SYNERGY SEATING SYSTEM CONTINUED****15.11. Joystick Bracket Receiver****REQUIRED.**

- ☐ Standard Fixed Receiver ..... Standard  
Part: FRMASMB16062
- ☐ Stealth Height Adjustable Receiver ..... \$214  
Tube Arm Part: ST-ARM250C-Q. HCPCS: K0108

**16. TRU-BALANCE FLEX SEATING SYSTEM**

**Optional.** Seat width and depth can be adjusted in 1" increments from 15"W-22"W and 16"D-20"D. Maximum seat size 22"W x 20"D. **Not available with power adjustable seat height options.**

**16.1. TRU-Balance Flex Frame Selection**

TRU-Balance Flex Seating System ..... No Charge  
Part: PTOINDV3853

**16.2. TRU-Balance Flex Seat Dimension - Width**

**REQUIRED with TRU-Balance Flex seating.** Maximum seat width is 22". Part number is SETWDTH plus number following # below.

|                  |                  |                  |                  |                  |                  |                  |                  |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| 15" W<br>(#1005) | 16" W<br>(#1006) | 17" W<br>(#1007) | 18" W<br>(#1008) | 19" W<br>(#1009) | 20" W<br>(#1010) | 21" W<br>(#1011) | 22" W<br>(#1012) |
|------------------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|

**16.3. TRU-Balance Flex Seat Dimension - Depth**

**REQUIRED with TRU-Balance Flex seating.** Maximum seat depth is 20". Part number is SETDPH plus number following # below.

|                  |                  |                  |                  |                  |                  |
|------------------|------------------|------------------|------------------|------------------|------------------|
| 15" D<br>(#1005) | 16" D<br>(#1006) | 17" D<br>(#1007) | 18" D<br>(#1008) | 19" D<br>(#1009) | 20" D<br>(#1010) |
|------------------|------------------|------------------|------------------|------------------|------------------|

**16.4. TRU-Balance Flex Seat to Floor**

**REQUIRED with TRU-Balance Flex seating.**

|                                              |                                                             |                                                              |
|----------------------------------------------|-------------------------------------------------------------|--------------------------------------------------------------|
| 16.5" at front and rear<br>Part: PTOINDV3459 | 17.5" at front with 2° posterior angle<br>Part: PTOINDV3460 | 18.25" at front with 4° posterior angle<br>Part: PTOINDV3805 |
|----------------------------------------------|-------------------------------------------------------------|--------------------------------------------------------------|

**16.5. TRU-Balance Flex Back Cane Selection**

Angle Adjustable Back Cane ..... No Charge  
Part: MEC140266

Flip Forward Back Cane ..... \$150  
Part: MEC140272

Not available with oxygen holder, or cane mounted and cantilever armrests.

2" Mesomorph (1" per side) Back Cane ..... \$195  
Part: MEC1804736

4" Mesomorph (2" per side) Back Cane ..... \$195  
Part: MEC1804739

2" Endomorph (1" per side) Back Cane ..... \$195  
Part: MEC141455

4" Endomorph (2" per side) Back Cane ..... \$195  
Part: MEC141456

Manual Recline Back ..... \$750  
Part: FRMASMB12910. HCPCS: E1226

Select Synergy Back and an overall back height, or to omit back. Back canes will be 22" H. Not available with flip back armrests, cantilever armrests, TRU-Comfort Back, Sport Back, rear accessories, or articulating foot platform. Pivot point is 2.5" above the seat pan.

**16.6. TRU-Balance Flex Back Cane Height Selection**

16" and 17" back canes are not available with a Synergy Back and cane mounted, flip up armrests.

14" High Canes  
Part: FRMASMB15832

17" High Canes  
Part: FRMASMB15835

20" High Canes  
Part: FRMASMB15838

23" High Canes  
Part: FRMASMB15841

15" High Canes  
Part: FRMASMB15833

18" High Canes  
Part: FRMASMB15836

21" High Canes  
Part: FRMASMB15839

24" High Canes  
Part: FRMASMB15842

16" High Canes  
Part: FRMASMB15834

19" High Canes  
Part: FRMASMB15837

22" High Canes  
Part: FRMASMB15840

25" High Canes  
Part: FRMASMB15843

**16.7. TRU-Balance Flex Back Options**

**REQUIRED with TRU-Balance Flex seating.** See the seat back availability chart for stock order back size ranges. Other size backs may be available with special order and extended lead time. Select a back type and overall back height for configured part number.

**1. Back Selection**

Synergy Back Cushion ..... \$275  
Part: See XRef

Sport Back Cushion, 20" Width or Less ..... \$375  
Part: See XRef

Sport Back Cushion, 22" Width ..... \$575  
Part: See XRef

TRU-Comfort Cushion, 20" Width or Less ..... \$675  
Part: See XRef. HCPCS: E2620

TRU-Comfort Cushion, 22" Width ..... \$775  
Part: See XRef. HCPCS: E2621

Sling Back ..... \$98  
Part: See XRef

Overall back height must match back cane height selection above.

Omit Back ..... No Charge  
Part: PTOINDV1287

|                |     | Back Width |     |     |     |     |     |     |     |   |   |   |   |   |   |   |   |
|----------------|-----|------------|-----|-----|-----|-----|-----|-----|-----|---|---|---|---|---|---|---|---|
|                |     | 15"        | 16" | 17" | 18" | 19" | 20" | 21" | 22" |   |   |   |   |   |   |   |   |
| Overall Height | 16" | -          | P   | T   | P   | T   | P   | T   | P   | T | P | T |   |   |   |   |   |
|                | 17" | S          | S   | P   | T   | S   | P   | T   | S   | P | T | S | P | T |   |   |   |
|                | 18" | S          | S   | P   | T   | S   | P   | T   | S   | P | T | S | P | T |   |   |   |
|                | 19" | S          | S   | P   | T   | P   | T   | S   | P   | T | S | P | T | S | P | T |   |
|                | 20" | S          | S   | P   | T   | P   | T   | S   | P   | T | S | P | T | S | P | T |   |
|                | 21" | S          | S   | P   | T   | S   | P   | T   | S   | P | T | S | P | T | S | P | T |
|                | 22" | S          | S   | P   | T   | S   | P   | T   | S   | P | T | S | P | T | S | P | T |
|                | 23" | S          | S   | P   | T   | S   | P   | T   | S   | P | T | S | P | T | S | P | T |
|                | 24" | S          | S   | P   | T   | S   | P   | T   | S   | P | T | S | P | T | S | P | T |
|                | 25" | S          | S   | P   | T   | S   | P   | T   | S   | P | T | S | P | T | S | P | T |

S

Synergy Back

P

Sport Back

T

TRU-Comfort Back

TRU-Balance Flex Seat Back Availability Chart

TRU-Balance Flex seating section continued on following pages.

**16. TRU-BALANCE FLEX SEATING SYSTEM CONTINUED****2. TRU-Balance Flex Overall Back Height**

Refer to the seat back availability chart for stock order back size ranges.

|                                        |                                        |                                        |                                        |                                        |
|----------------------------------------|----------------------------------------|----------------------------------------|----------------------------------------|----------------------------------------|
| <b>14" Height</b><br>Part: SETHGHT1027 | <b>17" Height</b><br>Part: SETHGHT1030 | <b>20" Height</b><br>Part: SETHGHT1033 | <b>23" Height</b><br>Part: SETHGHT1036 | <b>26" Height</b><br>Part: SETHGHT1039 |
| <b>15" Height</b><br>Part: SETHGHT1028 | <b>18" Height</b><br>Part: SETHGHT1031 | <b>21" Height</b><br>Part: SETHGHT1034 | <b>24" Height</b><br>Part: SETHGHT1037 | <b>27" Height</b><br>Part: SETHGHT1040 |
| <b>16" Height</b><br>Part: SETHGHT1029 | <b>19" Height</b><br>Part: SETHGHT1032 | <b>22" Height</b><br>Part: SETHGHT1035 | <b>25" Height</b><br>Part: SETHGHT1038 | <b>28" Height</b><br>Part: SETHGHT1041 |

**16.8. TRU-Balance Flex Armrests****REQUIRED with TRU-Balance Flex seating.** Left (L) and right (R) boxes may be checked to configure different size arms within one section arm type.**1a. TRU-Balance Flex Seat, Drop-In Armrests**

Includes full length, flat armpads.

|                                                                                                                                                                                                   |                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> Full Length, Drop-In Armrests ..... No Charge<br><input type="checkbox"/> <input type="checkbox"/> Left part: ARM1803253. Right part: ARM156742 | <input type="checkbox"/> <input type="checkbox"/> Desk Length, Drop-In Armrests ..... No Charge<br><input type="checkbox"/> <input type="checkbox"/> Left part: ARM1803254. Right part: ARM156743 |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**1b. Flip Back, Cane Mounted Armrests**

May interfere with swing-away laterals.

|                                                                                                                                                                                                          |                                                                                                                                                                                                          |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> Full Length, Flip Back Armrests (pair) ..... \$60<br><input type="checkbox"/> <input type="checkbox"/> Left part: FRMASMB7038. Right part: FRMASMB7039 | <input type="checkbox"/> <input type="checkbox"/> Desk Length, Flip Back Armrests (pair) ..... \$60<br><input type="checkbox"/> <input type="checkbox"/> Left part: FRMASMB7040. Right part: FRMASMB7041 |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**1c. Two-Post, Flip Back, Height Adjustable, Removable Armrests**

Not available with a depth less than 16". Required with flip forward back and laptrays.

|                                                                                                                                                                                    |                                                                                                                                                                                    |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> Full Length, Two-Post Armrests ..... \$150 ea.<br><input type="checkbox"/> <input type="checkbox"/> Part: See XRef. HCPCS: E0973 | <input type="checkbox"/> <input type="checkbox"/> Desk Length, Two-Post Armrests ..... \$150 ea.<br><input type="checkbox"/> <input type="checkbox"/> Part: See XRef. HCPCS: E0973 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**1d. Single Post, Height Adjustable, Removable Armrests**

Height adjustment from 11"-17" in half-inch increments.

|                                                                                                                                                                                                                       |                                                                                                                                                                                                                       |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> Full Length, Single Post Armrests ..... \$75 ea.<br><input type="checkbox"/> <input type="checkbox"/> Left part: FRMASMB8594. Right part: FRMASMB8595. HCPCS: E0973 | <input type="checkbox"/> <input type="checkbox"/> Desk Length, Single Post Armrests ..... \$75 ea.<br><input type="checkbox"/> <input type="checkbox"/> Left part: FRMASMB8592. Right part: FRMASMB8593. HCPCS: E0973 |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**16.9. TRU-Balance Flex Armpads****REQUIRED with armrest selection above except Drop-In Armrests.** Waterfall armpads may cause interference with the seat back. Armpads are configured as pairs.**1a. Straight Armpads**

|                                                                                  |                                                                                  |
|----------------------------------------------------------------------------------|----------------------------------------------------------------------------------|
| Straight, Full Length Left and Right ..... No Charge<br>Part: ARM134275          | Straight, Desk Length Left, Full Length Right ..... No Charge<br>Part: ARM134275 |
| Straight, Full Length Left, Desk Length Right ..... No Charge<br>Part: ARM134281 | Straight, Desk Length Left and Right ..... No Charge<br>Part: ARM134254          |

**1b. Waterfall Armpads**


|                                                                                   |                                                                                   |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
| Waterfall, Full Length Left and Right ..... No Charge<br>Part: ARM134275          | Waterfall, Desk Length Left, Full Length Right ..... No Charge<br>Part: ARM134275 |
| Waterfall, Full Length Left, Desk Length Right ..... No Charge<br>Part: ARM134281 | Waterfall, Desk Length Left and Right ..... No Charge<br>Part: ARM134254          |

**1c. Falcon Gel Armpads**

|                                                                                   |                                                                                   |
|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
| Falcon Gel, Full Length Left and Right ..... \$296.40<br>Part: ARM136438          | Falcon Gel, Desk Length Left, Full Length Right ..... \$274.22<br>Part: ARM136439 |
| Falcon Gel, Full Length Left, Desk Length Right ..... \$274.22<br>Part: ARM136442 | Falcon Gel, Desk Length Left and Right ..... \$252.40<br>Part: ARM136437          |

**17. HEADRESTS****Optional.** Select a headrest pad and bracket. **Not available with contoured seating.****17.1. Headrest Pads**

|                                                                                                                           |                                                                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> Stealth 8" Comfort Plus (CP450) ..... \$151.94<br>Part: POS148484. HCPCS: E0955                  | <input type="checkbox"/> Stealth Small Contoured (9.75"Wx3"H) ..... \$64.20<br>(STL-OTSM250)<br>Part: SETHEAD1052. HCPCS: E0955   |
| <input type="checkbox"/> Stealth 10" Comfort Plus (CP250) (A) ..... \$151.94<br>Part: SETHEAD1057. HCPCS: E0955           | <input type="checkbox"/> Stealth Large Contoured (11"Wx6"H) (B) ..... \$80.25<br>(STL-OTLG150)<br>Part: SETHEAD1053. HCPCS: E0955 |
| <input type="checkbox"/> Stealth 14" Comfort Plus (CP150) ..... \$151.94<br>Part: INDPART2327. HCPCS: E0955               |                                                                                                                                   |
| <input type="checkbox"/> Pediatric Comfort Plus (6"Wx4.5"H) ..... \$151.94<br>(CPN350)<br>Part: SETHEAD1069. HCPCS: E0955 |                                                                                                                                   |

**17.2. Headrest Brackets**

|                                                                                                                    |                                                                                                                                                                  |
|--------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Stealth's Worlds Best Removable Headrest Bracket (TWB485-1) (C) ..... \$220.50<br>Part: FRMASMB15999. HCPCS: E1028 | Stealth Multi-Axis Removable (D) ..... \$160.50<br>Headrest Bracket (HMO475P)<br>Part: HDWASMB1084. HCPCS: E1028<br>Not available with TRU-Balance Flex seating. |
| Pediatric Offset (HMO475P-17) (E) ..... \$160.50<br>Part: HDWASMB1166. HCPCS: E1028                                |                                                                                                                                                                  |

Headrest section continued on following page.

**17. HEADRESTS CONTINUED****17.3. Flip-down Headrest Bracket**

Stealth Flip-down Headrest Bracket (FDM350).....\$224.70  
 Part: INDPART2219  
 Not available with TRU-Balance Flex seating. The flip-down bracket is used with a headrest bracket above but will not be mounted to the chair due to permanent headrest post modification necessary to allow the headrest bracket to flip.

**18. LAP BELTS**

**REQUIRED except with select contoured seats with included lap belt.** Select one lap belt below.

**18.1a. Standard Lap Belts**

☐ 38" Lap Belt Part: ACCBELT1006 ☐ 50" Lap Belt Part: ACCBELT1000 ☐ 60" Lap Belt Part: ACCBELT1002 ☐ 70" Lap Belt Part: ACCBELT1003 ☐ 80" Lap Belt Part: ACCBELT1004 ☐ 90" Lap Belt Part: ACCBELT1001

**18.1b. Optional Lap Belts**

Not available with contoured or TRU-Balance Flex seating.

☐ i-Fit 60" Single Pull, Push Button, Medium (7.5") Padded, 1.5" Belt (IFPB5213) Part: ST-IFPB5213-Q .....\$86  
☐ i-Fit 60" Single Pull, Push Button, Extra Large Padded (11.5"), 2" Belt (IFPB5315) Part: ST-IFPB5315-Q .....\$86  
☐ i-Fit 60" Single Pull, Push Button, Large Padded (9.5"), 2" Belt (IFPB5314) Part: ST-IFPB5314-Q .....\$86  
☐ 60" Sure-Lok Lap Belt (Non-retractable).....\$120  
 Required with WC-19 Transit Option, Synergy Seat only.  
 Part: ACCBELT1016

**19. LATERAL EXTENSIONS**

**REQUIRED with Synergy Back and thoracic lateral selections.**

Lateral Extension for Synergy Back (pair).....\$26.25 ea.  
 16"OH back part: FRMASMB16274. 17"OH back part: FRMASMB16280. 20"OH back part: FRMASMB16298. 24"OH back part: FRMASMB16281

**20. SYNERGY POWER OPTION ELECTRONICS**

**REQUIRED for Synergy Seat with Power Legrests or Elevating Seat.** Select the desired control type for the seating system.

☐ One Actuator Control through NE+ Joystick .....\$1,895  
 Part: See XRef. HCPCS: E2310  
☐ One Actuator Control through Q-Logic 3 Electronics .....\$1,895  
 Part: See XRef. HCPCS: E2310  
☐ Two Actuator Control through NE+ Joystick .....\$2,375  
 Part: See XRef. HCPCS: E2311  
☐ Multiple Actuator Control through Q-Logic 3 Electronics .....\$2,895  
 Part: See XRef. HCPCS: E2311

**21. LEGRESTS**

**Optional.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference.

**21.1. Center Mount Foot Platform**

If power articulating foot platform is selected, be sure power legrest control has been selected in earlier section. Select legrest length in following section. **Not available with contoured seating.**

**1. Foot Platform Assembly**

Power Articulating Foot Platform .....\$2,995  
 Part: PTOINDV3754. HCPCS: E1012+K0040  
 Not available with TRU-Balance Flex seating.

Manually Adjustable Center Mount Foot Platform...\$250  
 12"-16" width part: FRMASMB17944. 14"-18" width part: FRMASMB17945. 16"-20" width part: FRMASMB17946.  
 HCPCS: K0040

**2. Lower Extension and Angle Adjustable Footplates**

8"-12" Extension, 5"x6" Footplates  
 Part: RIG152577  
 Not available with power articulating platform.

10"-16" Extension, 5"x6" Footplates  
 Part: RIG150206

13"-19" Extension, 5"x6" Footplates  
 Part: RIG151505

8"-12" Extension, 5"x8" Footplates  
 Part: RIG152578  
 Not available with power articulating platform.

10"-16" Extension, 5"x8" Footplates  
 Part: RIG150203

13"-19" Extension, 5"x8" Footplates  
 Part: RIG151506

10"-16" Extension, 5"x11" Footplates  
 Part: RIG150212

13"-19" Extension, 5"x11" Footplates  
 Part: RIG151508

| Center Mount Foot Platform<br>Compatibility Matrix | Synergy Seat |        |      |                                                                                                                                                                                                                                                    |
|----------------------------------------------------|--------------|--------|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                    | Low          | Medium | High |                                                                                                                                                                                                                                                    |
| 10"-16" Length Extension                           | !            | ✓      | ✓    | <input checked="" type="checkbox"/> No interference issues at this seat-to-floor height.<br><input type="checkbox"/> May have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability. |
| 13"-19" Length Extension                           | !            | !      | !    |                                                                                                                                                                                                                                                    |

**3. Center Mount Calf Support**

\*Calf pads are standard with a power articulating foot platform. Select calf pads covers for additional charge. Calf pad and panel covers include 1/2" of soft visco foam.

Small Calf Pads (pair) 6"Wx6"H .....\$150\*  
 Part: RIGASMB7120018

Small Calf Panel 18"Wx8"H .....\$150  
 Part: FRMASMB7110012

Large Calf Pads (pair) 8" taper to 6"Wx8"H .....\$150\*  
 Part: RIGASMB7120019

Large Calf Panel 24"Wx8"H .....\$200  
 Part: RIG123824  
 Only available on widths of 20" or greater.

Legrests section continued on following page.

**21. LEGRESTS CONTINUED****4. Center Mount Modular Footplate Extension Components**

Optional. Assembly required. May require seating configuration adjustments to avoid front caster interference.

|                                                                                                                                                             |                                                                                                                             |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Width Extension for 8"D Footplates .....\$50<br>Part: RIG152517                                                                                             | 4" Outside Edge (pair) .....\$50<br>Part: RIG152510<br>Only available with 6"D footplates and 8"D footplate extension.      |
| Width Extension for 11"D Footplates .....\$50<br>Part: RIG152518                                                                                            | 6" Outside Edge (pair) .....\$50<br>Part: RIG152516<br>Available with 8"D and 11"D footplates and 11"D footplate extension. |
| 2" Bell-Shaped Length Extension .....\$50<br>Part: RIG152519<br>Not compatible with width extensions. Can be used with 6" outside edge with 11" footplates. |                                                                                                                             |

**5. Center Mount Additional Accessories**

|                                                                                                   |                                                               |
|---------------------------------------------------------------------------------------------------|---------------------------------------------------------------|
| Reversed Dartex Calf Pad Cover (pair) ..... \$150<br>Small part: RIG133717. Large part: RIG133719 | Wheel Bumpers (pair)..... \$12.50 ea.<br>Part: WHLASMB1578 x2 |
| Reversed Dartex Panel Cover (pair) ..... \$150<br>Small part: RIG134673. Large part: RIG134744    | Footplate Connector ..... \$15<br>Part: RIGASMB120023         |

**6. Stealth Gel Footplate Covers**

|                                                                                     |                                                                                                          |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------|
| 5"x8" Tapered Gel Footplate Covers..... \$205<br>Part: ST-FPG5008T-Q. HCPCS: K0108  | 6.5"x8" Tapered Gel Footplate Covers for Width Extension..... \$225<br>Part: ST-FPG6508T-Q. HCPCS: K0108 |
| 5"x11" Tapered Gel Footplate Covers..... \$225<br>Part: ST-FPG5011T-Q. HCPCS: K0108 | 7"x11" Tapered Gel Footplate Cover for Width Extension ..... \$250<br>Part: ST-FPG7011T-Q. HCPCS: K0108  |

**21.2. 70° Power Articulating, Swing Away Legrests**

**Not available with contoured or TRU-Balance Flex seating.** 4" of articulating range. Be sure power legrest control has been selected in previous section. Select an upper extension type and footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Articulating Legrest Upper Assembly**

|                                                                                                                                                                                                                                                                                 |                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 70° Articulating, Swing-Away.....\$1,197.50 ea.<br>Legrest with 14.5"-18" Extension<br>Pair Part: FRMASMB9073. HCPCS for pair: E1010+K0040+E0951.<br>Left Part: FRMASMB9074. Right Part: FRMASMB9075 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 70° Articulating, Swing-Away.....\$1,197.50 ea.<br>Legrest with 11"-13" Extension<br>Pair Part: FRMASMB11444. HCPCS for pair: E1010+K0040+E0951.<br>Left Part: FRMASMB11366. Right Part: FRMASMB11367. |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

| 70° Power Articulating, Swing-Away<br>Legrest Compatibility Matrix | Synergy Seat |        |      |                                                                                                                                                                                                                                                                                                                 |
|--------------------------------------------------------------------|--------------|--------|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                    | Low          | Medium | High |                                                                                                                                                                                                                                                                                                                 |
| Angle Adjustable, 11" - 13"                                        | ✓            | ✓      | ✓    | <input checked="" type="checkbox"/> This rigging configuration causes no interference issues at this seat-to-floor height.<br><input type="checkbox"/> This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability. |
| Angle Adjustable, 14.5" - 18"                                      | !            | !      | !    |                                                                                                                                                                                                                                                                                                                 |

**2. Angle Adjustable Footplate(s) with Heel Loop and Wheel Bumper**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

|                                                                                                                                                            |                                                                                                                                                                         |                                                                                                                                                                          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D<br>Pair part: FRMASMB7448. Left part: RIG133543. Right part: RIG133566. | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx8"D<br>Pair part: FRMASMB7450. Left part: RIG133545. Right part: RIG133568.              | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D <sup>1</sup><br>Pair part: FRMASMB7452. Left part: RIG133548. Right part: RIG133570. |
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D<br>Pair part: FRMASMB7449. Left part: RIG133544. Right part: RIG133567. | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx8"D <sup>1</sup><br>Pair part: FRMASMB7451. Left part: RIG133547. Right part: RIG133569. | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D <sup>2</sup><br>Pair part: FRMASMB7453. Left part: RIG133549. Right part: RIG133571. |

**3. Power Articulating, Swing-Away Legrest Accessories**

|                                                           |                                                                                                                                                                                                         |
|-----------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Calf Pad (pair) 6"Wx6"H..... No Charge<br>Part: MEC153735 | Falcon Gel Padded Knee Adductor Buttons (KBE)(pair) .....\$125<br>Part: INDPART2242. HCPCS: E0953^<br>^Please verify the appropriate billing code to be used for this item with all third-party payors. |
|-----------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**21.3. Swing-Away Legrests**

Select a legrest upper, a lower extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes.

**1a. 60°-90° Swing-Away Legrest Uppers**

Adjustable in 5° increments.

|                                                                                                                                                                                                                                   |                                                                                                                                                                                                                                           |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 60°-90° Swing-Away Legrest..... \$81 ea.<br>Pair part: RIG120062. Left: FRMASMB17558.<br>Right: FRMASMB17559<br>Not available with 11"-13" extensions. | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 60°-90° Swing-Away Legrest, Short ..... \$81 ea.<br>Pair part: RIG136481. Left: RIG132607.<br>Right: RIG132609.<br>For use only with 11"-13" lower extensions. |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**1b. 70° Swing-Away Legrest Uppers**

Only available with occupied transit WC-19 selections. Sold only in pairs. Select lower extension and footplates as pairs.

|                                                                                                                       |                                                                                                                                  |
|-----------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| 70° Swing-Away Legrest (pair).....\$162<br>Pair part: FRMASMB7841<br>Not available with 11"-13" extensions.           | 70° Tapered, Swing-Away Legrest, Short (pair).....\$162<br>Pair part: FRMASMB8052<br>For use only with 11"-13" lower extensions. |
| 70° Tapered, Swing-Away Legrest (pair) .....\$162<br>Pair part: FRMASMB7843<br>Not available with 11"-13" extensions. |                                                                                                                                  |

Legrests section continued on following page.

**21. LEGRESTS CONTINUED**

| Swing-Away Legrest Compatibility Matrix | Contoured/Synergy Seat |        |      |                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------------------|------------------------|--------|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                         | Low                    | Medium | High |                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Extruded, 13" - 16.5"                   | !                      | ✓      | ✓    | <input checked="" type="checkbox"/> This rigging configuration causes no interference issues at this seat-to-floor height.<br><input type="checkbox"/> This rigging configuration may have less than 1.75" clearance when fully extended and may interfere with the ground or reduce obstacle climbing ability.<br><input checked="" type="checkbox"/> This rigging configuration does not work with this seat-to-floor height. |
| Extruded, 15.5" - 19.5"                 | !                      | !      | !    |                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Angle Adjustable, 13" - 16.5"           | !                      | ✓      | ✓    |                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Angle Adjustable, 15.5" - 19.5"         | !                      | !      | !    |                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Multi-Axis, 15.5" - 19.5"               | !                      | !      | !    |                                                                                                                                                                                                                                                                                                                                                                                                                                 |

**2. Lower Extension**

|                                                                                                                                     |                                                                                                                                             |                                                                                                                                       |                                                                                                                 |
|-------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 11"-13" Length<br>Part: See XRef<br>N/A with Multi-axis. | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 11"-16.5" Length<br>Part: See XRef<br>N/A with Angle Adjustable. | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 13"-16.5" Length<br>Part: See XRef<br>N/A with Multi-Axis. | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 15.5"-19.5" Length<br>Part: See XRef |
|-------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|

**3a. Extruded Footplate(s)**

|                                                                                                                                                                                           |                                                                                                                                                                                                                                                   |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 7.25"Wx6"D Extruded Footplate..... Standard<br>Part pair: FRMASMB7847. Left: FRMASMB7861<br>Right: FRMASMB7862 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8.5"Wx6"D Extruded Footplate ..... Standard<br>Part pair: FRMASMB7848. Left: FRMASMB7863<br>Right: FRMASMB7864<br>May cause interference on seat widths less than 18." |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**3b. Angle Adjustable Footplate(s) with Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

|                                                                                                                                                                                       |                                                                                                                                                                                                     |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D ..... \$85 ea.<br>Pair part: FRMASMB10564. Left: RIG133537. Right: RIG133560<br>HCPCS: K0040+E0951 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx8"D <sup>1</sup> ..... \$85 ea.<br>Pair part: FRMASMB10567. Left: RIG133540. Right: RIG133563<br>HCPCS: K0040+E0951  |
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D..... \$85 ea.<br>Pair part: FRMASMB10565. Left: RIG133538. Right: RIG133561<br>HCPCS: K0040+E0951  | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D <sup>1</sup> ..... \$85 ea.<br>Pair part: FRMASMB10568. Left: RIG133541. Right: RIG133564<br>HCPCS: K0040+E0951 |
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx8"D ..... \$85 ea.<br>Pair part: FRMASMB10566. Left: RIG133539. Right: RIG133562<br>HCPCS: K0040+E0951 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D <sup>2</sup> ..... \$85 ea.<br>Pair part: FRMASMB10569. Left: RIG133542. Right: RIG133565<br>HCPCS: K0040+E0951 |

**3c. Multi-Axis Footplate(s) with Heel Loops**

<sup>1</sup>The 6"x11" footplates may cause interference with seat widths less than 18". <sup>2</sup>The 8"x11" footplates may cause interference with seat widths less than 20".

|                                                                                                                                                                                          |                                                                                                                                                                                                        |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D ..... \$247.50 ea.<br>Pair part: FRMASMB9456. Left: RIG133531. Right: RIG133552<br>HCPCS: K0040+E0951 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx8"D <sup>1</sup> ..... \$247.50 ea.<br>Pair part: FRMASMB9459. Left: RIG133534. Right: RIG133555<br>HCPCS: K0040+E0951  |
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D..... \$247.50 ea.<br>Pair part: FRMASMB9457. Left: RIG133532. Right: RIG133553<br>HCPCS: K0040+E0951  | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 6"Wx11"D <sup>1</sup> ..... \$247.50 ea.<br>Pair part: FRMASMB9460. Left: RIG133535. Right: RIG133558<br>HCPCS: K0040+E0951 |
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx8"D ..... \$247.50 ea.<br>Pair part: FRMASMB9458. Left: RIG133533. Right: RIG133554<br>HCPCS: K0040+E0951 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 8"Wx11"D <sup>2</sup> ..... \$247.50 ea.<br>Pair part: FRMASMB9461. Left: RIG133536. Right: RIG133559<br>HCPCS: K0040+E0951 |

**4. Swing-Away Legrest Accessories**

|                                                                                                                           |                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|
| Heel Loops (pair)..... \$50<br>Part: FRMASMB7873. HCPCS: E0951<br>Standard on angle adjustable and multi-axis footplates. | Wheel Bumper..... \$12.50 ea.<br>Part: See XRef                                   |
|                                                                                                                           | Gel Neoprene Footrest Sleeve (pair) ..... \$80<br>Part: INDPART2243. HCPCS: K0108 |

**5. Amputee Supports with Swing-Away Upper**

Enter a quantity for the amputee supports and covers. A 60°-90° legrest upper must be selected with the amputee support. The support will be configured with the legrest upper side that does not have a lower extension and footplate selected with it.

|                                                                                              |                                                                                             |
|----------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
| Therafin 31515 Multi-Axis Amputee Support ..... \$230 ea.<br>Part: INDPART2822. HCPCS: E1020 | Therafin 32387 Amputee Support Gel Cover ..... \$155 ea.<br>Part: INDPART2823. HCPCS: K0108 |
|----------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|

**21.4. Pediatric, Swing-Away Legrests**

Not available with TRU-Balance Flex seating. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Pediatric Legrests Upper**

|                                                                                                                                                                                                                           |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 90°, Pediatric, Swing-Away Legrest (7"-10" Length Extension) ..... \$81 ea.<br>Pair part: FRMASMB8808. Left: FRMASMB8809<br>Right: FRMASMB8810 |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**2. Angle Adjustable Footplate(s)**

|                                                                                                                                                                                                 |                                                                                                                                                                                                 |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 4"Wx6"D Angle Adjustable ..... \$85 ea.<br>Pair part: FRMASMB8812. Left: RIG133537. Right: RIG133560<br>HCPCS: K0040 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> 5"Wx6"D Angle Adjustable ..... \$85 ea.<br>Pair part: FRMASMB8811. Left: RIG133538. Right: RIG133561<br>HCPCS: K0040 |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**3. Accessories**

|                                                                |                                                                                                                                                                          |
|----------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Heel Loops (pair)..... \$50<br>Part: FRMASMB7873. HCPCS: E0951 | <input type="checkbox"/> <input type="checkbox"/> <input type="checkbox"/> Wheel Bumper(s) ..... \$12.50 ea.<br>Pair part: WHL169508. Left: WHL169510. Right: WHL169511. |
|----------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

Legrests section continued on following page.

**21. LEGRESTS CONTINUED****21.5. Elevating Legrests**

Select a legrest upper, a legrest extension, and type of footplate. Pair (P), left (L), and right (R) configurations are selected by checking the appropriate boxes. Select legrest length in following section.

**1. Elevating Legrests Upper**

☐ ☐ ☐ Elevating Legrest Assembly ..... \$162.50 ea.  
 Pair part: FRMASMB6483. Left: RIG138012  
 Right: RIG138013. HCPCS: E0990  
 Select lower extension and footplates

**2. Lower Extension**

☐ ☐ ☐ 11"-14" Length ☐ ☐ ☐ 14"-17" Length  
 Pair: FRMASMB4755 Pair: FRMASMB4756  
 Left: RIG138026 Left: RIG138033  
 Right: RIG138030 Right: RIG138035  
 For reinforced only For reinforced only

☐ ☐ ☐ 12"-14.5" Length ☐ ☐ ☐ 14.5"-18" Length  
 Part: See XRef Part: See XRef  
 For angle adjustable only For angle adjustable only

**3a. Reinforced Cast Footplate(s)**

☐ ☐ ☐ 7.25"Wx6"D Reinforced Footplate..... Standard ☐ ☐ ☐ 8.5"Wx6"D Reinforced Footplate..... Standard  
 Pair part: FRMASMB4748. Left: RIG138055. Right: RIG138056 Pair part: FRMASMB4747. Left: RIG138060. Right: RIG138062  
 May cause interference on seat widths less than 18."

**3b. Angle Adjustable Footplate(s)**

\*The 6"x11" footplates may cause interference with seat widths less than 18". \*The 8"x11" footplates may cause interference with seat widths less than 20". 4"x6" and 5"x6" footplates are only available in pairs.

☐ ☐ ☐ 4"Wx6"D ..... \$85 ea.  
 Pair parts: FRMASMB8812 x2.  
 HCPCS: K0040

☐ ☐ ☐ 5"Wx6"D..... \$85 ea.  
 Pair parts: FRMASMB8811 x2.  
 HCPCS: K0040

☐ ☐ ☐ 5"Wx8"D Footplates with Heel Loops ..... \$85 ea.  
 Pair part: FRMASMB4735. Left: RIG137791. Right: RIG137792  
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx8"D<sup>1</sup> Footplates with Heel Loops ..... \$85 ea.  
 Pair part: FRMASMB4738. Left: RIG137793. Right: RIG137795  
 HCPCS: K0040+E0951

☐ ☐ ☐ 6"Wx11"D<sup>1</sup> Footplates with Heel Loops ..... \$85 ea.  
 Pair part: FRMASMB4741. Left: RIG137796. Right: RIG137797  
 HCPCS: K0040+E0951

☐ ☐ ☐ 8"Wx11"D<sup>2</sup> Footplates with Heel Loops..... \$85 ea.  
 Pair part: FRMASMB4744. Left: RIG137798. Right: RIG137799  
 HCPCS: K0040+E0951

**4. Elevating Legrest Accessories**

Heel Loops (pair)..... \$50  
 Part: FRMASMB7873. HCPCS: E0951  
 Standard on angle adjustable footplates.

Wheel Bumper..... \$12.50 ea.  
 Part: See XRef  
 Available on angle adjustable footplates.

Falcon Gel Padded Knee Adductor Buttons (KBE) (pair) .....\$125  
 Part: INDPART2242. HCPCS: E0953^  
 ^Please verify the appropriate billing code to be used for this item with all third-party payors.

**22. SET LEGREST LENGTH**

**REQUIRED with matching left and right legrest configuration.** All efforts will be made to accommodate legrest length selection, but minor changes may need to be made to meet shipping guidelines. Some adjustments may be needed upon receipt of the unit. Reference the section's compatibility matrix for legrest interference. Non-matching or individual legrest configurations lengths are not set by the manufacturer.

|                   |                   |                    |                    |                    |                    |                    |                    |                    |                    |                    |
|-------------------|-------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|--------------------|
| 7" L<br>(#3470)   | 8.5" L<br>(#3473) | 10" L<br>(#3476)   | 11.5" L<br>(#3479) | 13" L<br>(#3482)   | 14.5" L<br>(#3485) | 16" L<br>(#3488)   | 17.5" L<br>(#3491) | 19" L<br>(#3494)   | 20.5" L<br>(#3497) | 22" L<br>(#3669)   |
| 7.5" L<br>(#3471) | 9" L<br>(#3474)   | 10.5" L<br>(#3477) | 12" L<br>(#3480)   | 13.5" L<br>(#3483) | 15" L<br>(#3486)   | 16.5" L<br>(#3489) | 18" L<br>(#3492)   | 19.5" L<br>(#3495) | 21" L<br>(#3498)   | 22.5" L<br>(#3670) |
| 8" L<br>(#3472)   | 9.5" L<br>(#3475) | 11" L<br>(#3478)   | 12.5" L<br>(#3481) | 14" L<br>(#3484)   | 15.5" L<br>(#3487) | 17" L<br>(#3490)   | 18.5" L<br>(#3493) | 20" L<br>(#3496)   | 21.5" L<br>(#3600) |                    |

**23. SEAT CUSHION**


**Optional.** Specific cushion/component sizes available, part numbers, and special order cushions can be found on the respective complete [Stealth Cushion Order Forms](#). This section allows a quick selection of a single cushion. The XXYY portion of the part number directly relates to the dimensions of the cushion, where XX represents the width and YY represents the depth.

**23.1. Seat Cushion Model Selection****Stealth Simplicity G Cushion**

General use cushion with Stealth Tek 2, 2-way stretch, breathable cover.  
Part series: *ST-SIMGXXYY*

Simplicity G, 10"W to 21" W, E2601 ..... \$110


Simplicity G, 22"W, E2602 ..... \$165

**Stealth Spectrum Foam SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: *ST-SPCFSPXXYY*

Spectrum Foam SPP, 16"W to 20" W, E2607 ..... \$460


Spectrum Foam SPP, 22"W, E2608 ..... \$530

**Stealth Solution SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: *ST-SOLSPPXXYY*

Solution SPP, 10"W to 21" W, E2607 ..... \$375


Solution SPP, 22"W, E2608 ..... \$430

**Stealth Glacial SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-GLSPXXYY*

Glacial SP, 10"W to 21"W, E2603 ..... \$330

Glacial SP, 22"W, E2604 ..... \$440

**Stealth Zen SP Cushion**

Skin protection with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-ZENSPXXYY*

Zen SP, 10"W to 21" W, E2603 ..... \$315


Zen SP, 22"W, E2604 ..... \$425

**Stealth TRU-Comfort 2 SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-TC2SPPXXYY*

TRU-Comfort 2 SPP, 16"W to 20" W, E2607 ..... \$400


TRU-Comfort 2 SPP, 22"W, E2608 ..... \$500

**Stealth Spectrum Gel SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover. Includes urethane liner.  
Part series: *ST-SPCGSPPXXYY*

Spectrum Gel SPP, 14"W to 20" W, E2607 ..... \$460

Spectrum Gel SPP, 22"W, E2608 ..... \$530

**Stealth Essence SPP Cushion**

Skin protection and positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-ESSSPPXXYY*

Essence SPP, No Wedge, 16"-20" Wide, E2607 ..... \$400

Essence SPP, No Wedge, 22" Wide, E2608 ..... \$500

Essence SPP, Small Wedge, 16"-20" Wide, E2607 .... \$435

Essence SPP, Small Wedge, 22" Wide, E2608 ..... \$535

Essence SPP, Medium Wedge, 16"-20" Wide, E2607 . \$465

Essence SPP, Medium Wedge, 22" Wide, E2608 ..... \$565

Essence SPP, Large Wedge, 16"-20" Wide, E2607 .... \$495

Essence SPP, Large Wedge, 22" Wide, E2608 ..... \$595

**Stealth Premiere P Cushion**

Positioning with Coolcore, 4-way stretch, breathable cover.  
Part series: *ST-PREMPXXYY*

Premiere P, 16"W to 20" W, E2605 ..... \$325

Premiere P, 22" W, E2606 ..... \$475

**23.2. Seat Cushion Dimensions**

**REQUIRED** with cushion model selection. See respective [cushion order form](#) for available sizes. Select to match or pick non-matching dimensions.

Match Seat Width and Depth

**Cushion Width**

|       |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| 10" W | 14" W | 16" W | 18" W | 20" W | 22" W |
| 12" W | 15" W | 17" W | 19" W | 21" W |       |

**Cushion Depth**

|       |       |       |       |       |       |
|-------|-------|-------|-------|-------|-------|
| 10" D | 14" D | 16" D | 18" D | 20" D | 22" D |
| 12" D | 15" D | 17" D | 19" D | 21" D |       |

**23.3. Seat Cushion Substitute Fluid Proof Liners**

Fluid Proof Urethane Liner ..... \$50-\$70

Part: URL1709070-XXYY

Not available with Simplicity, Essence, and TRU-Comfort 2.

Fluid Proof Polyurethane Liner ..... \$55-\$75

Essence/TC2 part: PUL1709080-XXYY. Others part: PUL1709030-XXYY

Not available with Simplicity.

**23.4. Seat Cushion Substitute Cover**

The standard cover may be substituted with a cover by checking the option below. The ""prefix" of the part number will be the part prefix of the cushion model selected above. For example, SPCGSPP-NBLK1816 will be an additional coolcore cover for a Spectrum Gel cushion.

Stealth Tek 2 Cover ..... No Charge

Part: \*prefix-TBLKXXYY

Not available with Simplicity, TRU-Comfort 2, Essence, or Premiere.

Mesh Cover ..... No Charge

Part: \*prefix-MBLKXXYY

Not available with Simplicity or TRU-Comfort 2.

**24. ACCESSORIES**

**Optional.** Multiple accessories can be mounted along the various tracks but some interference may occur. Accessories may need to be installed by the provider.

Contoured Seat Cup Holder ..... \$20

Part: ACCASMB2287. Part with 12" lift: ACCASMB2288

Synergy/Flex Seat Cup Holder ..... No Charge

Part: ACCASMB2634

Not available on TRU-Balance Flex seating without two-post armrests.

Oxygen Holder for Contoured Seat ..... \$118

Part: ACCASMB2292. HCPCS: E2208

Oxygen Holder for Synergy Seat ..... \$150

Not available with Synergy Manual Recline.

Part: FRMASMB7030. HCPCS: E2208

Oxygen Holder for Flex Seat ..... \$200

Part: ACC125006. HCPCS: E2208

Not available with Synergy Manual Recline.

Walker Holder ..... \$103

Contoured seat part: ACCASMB1912. Synergy/Flex part: FRMASMB7031

Cane/Crutch Holder ..... \$110

Contoured seat part: ACCASMB2294. HCPCS: E2207

Synergy part: FRMASMB4822. Flex part: ACC152515. HCPCS: E2207

Cell Phone Holder for 12" Lift Captain Seat ..... \$70

Part: ACC1611300

Accessories section continued on following page.

**24. ACCESSORIES CONTINUED**

|                                                                                                     |                                                                                                                        |
|-----------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|
| <b>Rear Basket</b> .....\$118<br>Contoured seat part: ACCASMB1914<br>Synergy part: FRMASMB4821      | <b>Quantum Backpack</b> ..... \$150<br>Part: ACC1905853                                                                |
| <b>Flag Holder</b> .....\$25<br>Contoured seat part: ACCASMB1930<br>Synergy/Flex part: ACCASMB2415. | <b>Grip Lap Board (20005)</b> .....\$38.95<br>Part: ACC152992                                                          |
| <b>Saddle Bag</b> .....\$35<br>Part: DWR1020D010                                                    | <b>Grip Drink Holder (20085)</b> .....\$15.50<br>Part: ACC152993                                                       |
| <b>Stealth Cell Phone/Media Device Holder</b> .....\$195<br>Part: ACC7120006                        | <b>Accessory Bar</b> ..... \$100<br>12"W-16"W part: MEC145955. 16"W-20"W part: MEC143975.<br>20"W-22"W part: MEC146018 |
| <b>Stealth Large Tablet Holder</b> ..... \$375<br>Part: ACC7120005                                  | <b>Large, Black, Weather Cover</b> ..... \$202<br>Part: ACCCOVR1021                                                    |

**25. BRANCH PATCHES**

Optional.

|                                                      |                                                        |                                                         |
|------------------------------------------------------|--------------------------------------------------------|---------------------------------------------------------|
| <b>Air Force Patch</b> ..... \$20<br>Part: ACC140240 | <b>Coast Guard Patch</b> ..... \$20<br>Part: ACC140239 | <b>Navy Patch</b> ..... \$20<br>Part: ACC140244         |
| <b>Army Patch</b> ..... \$20<br>Part: ACC140241      | <b>Marines Patch</b> ..... \$20<br>Part: ACC140243     | <b>Purple Heart Patch</b> ..... \$20<br>Part: ACC140238 |

**26. PROGRAMMER**

Optional.

|                                                                                           |                                                                                                         |
|-------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <b>Q-Logic EX/NE Programmer</b> ..... \$545<br>Part: ELEASMB6976 | <input type="checkbox"/> <b>Q-Logic PC Programming Disk with Cable</b> ..... \$750<br>Part: ELEASMB5216 |
|-------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|

**27. ALTERNATIVE CONTROLS AND ELECTRONICS**

Optional. Only available with the Synergy Seating system. Select a type of specialty controls below or to refer to the specialty controls order form.

**27.1. Enhanced Display**

**Q-Logic 3 Enhanced Display Kit**.....\$1,720  
Part: CTL167726  
Select a color option below. Includes enhanced display, flexible gooseneck and mounting bracket. Features: 3.5" LCD display; up to 5 profiles to control drive, seat, and auxiliary functions; Bluetooth connectivity for mobile devices and mouse emulation; integrated single switch scanner; IR controls.

**Q-Logic 3 Enhanced Display Shroud**

Required with Q-Logic 3 Enhanced Display selection above. Reference colors shown on page 2.

|                                                                      |                                                                                |
|----------------------------------------------------------------------|--------------------------------------------------------------------------------|
| <b>Black, Soft Touch</b> .....Standard<br>Part: DGN137728            | <b>Match Joystick Shroud to Base Color</b> ..... \$50<br>Part: See color below |
| <b>Non-matching Joystick Shroud</b> ..... \$50<br>Select color below |                                                                                |



|                                              |                                                |                                                 |                                               |                                                      |
|----------------------------------------------|------------------------------------------------|-------------------------------------------------|-----------------------------------------------|------------------------------------------------------|
| <b>Red Corvette</b><br>Part: DGN137729-PRD1  | <b>White Satin</b><br>Part: DGN137729-PWH2     | <b>Evolution Orange</b><br>Part: DGN137729-POR7 | <b>Lemon Crush</b><br>Part: DGN137729-PYL4    | <b>Military Digital Camo</b><br>Part: DGN137729-WCM7 |
| <b>Cherry Bomb</b><br>Part: DGN137729-PRD4   | <b>Deep Purple</b><br>Part: DGN137729-PPR5     | <b>Midnight Blue</b><br>Part: DGN137729-PBU10   | <b>Green Machine</b><br>Part: DGN137729-PGN11 | <b>Faded Glory</b><br>Part: DGN137729-WFL2           |
| <b>Back in Black</b><br>Part: DGN137729-PBK1 | <b>Raspberry Beret</b><br>Part: DGN137729-PPK4 | <b>Electric Blue</b><br>Part: DGN137729-PBU9    | <b>Military Green</b><br>Part: DGN137729-PGN8 | <b>Carbon Fiber</b><br>Part: DGN137729-WCB1          |
| <b>Night Train</b><br>Part: DGN137729-PBK3   | <b>Pretty in Pink</b><br>Part: DGN137729-PPK3  | <b>Blue Jean</b><br>Part: DGN137729-PBU11       | <b>Sugar Plum</b><br>Part: DGN137729-PPR7     | <b>To Be Determined</b><br>Follow-up required.       |
| <b>Grey Street</b><br>Part: DGN137729-PSL3   |                                                |                                                 |                                               |                                                      |

**27.2. SCIM Module**

Optional.

|                                                                                                                                                                                                                                                                                                                                                                                                          |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <b>SCIM Module for Q-Logic 3</b> .....\$525<br>Part: CTL168055<br>Allows use of joystick display for specialty devices. Available with Q-Logic 3 EX Joystick selection above. Up to 5 profiles to control drive, seat, and auxiliary functions, excluding sip-n-puff and single switch functions. A switch is required to operate the input device attached to the SCIM module. |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**27.3. Q-Logic 3 SAJ**

|                                                                                                                                                            |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <input type="checkbox"/> <b>Q-Logic 3 EX SAJ</b> .....\$700<br>Part: CTL143980<br>Must be used in conjunction with Enhanced Display or Q-Logic 3 Joystick. |
|------------------------------------------------------------------------------------------------------------------------------------------------------------|



Alternative controls and electronics section continued on following page.

**27. ALTERNATIVE CONTROLS AND ELECTRONICS CONTINUED****27.4. Head Array**

Adult, i-Drive Tri-Array Series, 3-Switch ..... \$4,414  
Head Array (IDHBT200-1)

Part: ST-IDHBT200-1-Q. HCPCS: E2330

Features bluetooth, head support, 3 non-proportional proximity sensors for drive control, an egg switch for reset mode change, and i-Drive smart box.

Adult, i-Drive Pro Series Head Array ..... \$4,461  
(IDHBT300-1)

Part: ST-IDHBT300-1-Q. HCPCS: E2330

Features bluetooth, head support with occipital pad and 2 swing-away spot pads for left and right side, 3 non-proportional proximity sensors, an egg switch for reset mode change, and i-Drive smart box.



Adult Pro-Series

**27.5. Sip and Puff Drive****1. Sip and Puff Module**

Sip and Puff Module for Q-Logic 3 ..... \$1,895

Part: CTL168056. E2325

Includes module and mounting hardware.

**2. Sip and Puff Straw and Mount**

Stealth Sip and Puff for Comfort Series, ..... \$294.25  
Left Mount (SP625-CL)

Part: ST-SP625-CL-Q. E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Comfort Series, ..... \$294.25  
Right Mount (SP622-CR)

Part: ST-SP625-CR-Q. E2326

Includes straw and mounting hardware. Use with Stealth Comfort Plus headrests.

Stealth Sip and Puff for Stealth i-Drive, ..... \$294.25  
Head Array (SP625-U)

Part: ST-SP625-U-Q. E2326

Includes straw and mounting hardware. Use with Stealth i-Drive Head Arrays.



Sip &amp; Puff Module



Stealth Sip &amp; Puff

**27.6. Egg Switches**

Red Egg Switch (ESRED) ..... \$80.85

Part: SWTMCR01077

Blue Egg Switch (ESBLU) ..... \$80.85

Part: SWTMCR01079

Black Egg Switch (ESBLK) ..... \$80.85

Part: SWTMCR01081

Green Egg Switch (ESGRN) ..... \$80.85

Part: SWTMCR01078

Yellow Egg Switch (ESYEL) ..... \$80.85

Part: SWTMCR01080



Egg Switch

**27.7. Switch Adapter**

2 to 1 Switch Adapter ..... \$13.39

Part: DWR1361H037

Connects 2 single switches to a stereo jack to act as a dual switch.

**27.8. Environmental Controls**

Q-Logic 3 ECU Module Kit ..... \$2,400

Part: CTL168057

**28. POSITIONING COMPONENTS**

Optional. Not available with contoured seating. More options are available on the [Positioning Components Order Form](#).

**28.1. Stealth Thoracic Laterals**

Swing-away laterals may interfere with TRU-Balance 3 armrests and 2-post, flip up armrests. These laterals are not available with the back cane option.

**1. Bracket**

Swing-away, Rail Mounts, pair (TWBSTL) ..... \$472.50

Part: INDPART2896. HCPCS: E1028

**2. Pads**

Curved Pads, 3"x5", pair (TWBL-3X5) ..... \$210

Part: INDPART2663. HCPCS: E0956

Curved Pads, 4"x6", pair (TWBL-4X6) ..... \$210

Part: INDPART2664. HCPCS: E0956

**3. Installation**

Install Laterals  
Part: PTOINDV3499

Do Not Install  
Part: PTOINDV3500

**28.2. Stealth UniLink Pelvic/Thigh Guide****1. UniLink Pelvic/Thigh Hardware**

UniLink, 5" Height Adjustable, ..... \$245 ea.  
Removable Mounting Hardware

Part: ST-UTTP5-Q

UniLink, 7" Height Adjustable, ..... \$255 ea.  
Removable Mounting Hardware

Part: ST-UTTP7-Q

**2. UniLink Pelvic/Thigh Pads**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

UniLink, 4" x 4" Pad ..... \$95 ea.

Part: ST-UTP544056100-Q

UniLink, 4" x 12" Pad ..... \$125 ea.

Part: ST-UTP54D056100-Q

UniLink, 4" x 8" Pad ..... \$105 ea.

Part: ST-UTP548056100-Q



UniLink Pad &amp; Hardware

Positioning components section continued on following page.

**28. POSITIONING COMPONENTS CONTINUED****3. UniLink Pelvic/Thigh Additional Links**

|                                                               |                                                           |
|---------------------------------------------------------------|-----------------------------------------------------------|
| UniLink, 1.5" Link ..... \$75 ea.<br>Part: ST-UIL001150A000-Q | UniLink, 2" Link..... \$75 ea.<br>Part: ST-UIL001200A00-Q |
|---------------------------------------------------------------|-----------------------------------------------------------|

**4. UniLink Pelvic/Thigh Guide Installation**

|                                      |                                       |                                     |
|--------------------------------------|---------------------------------------|-------------------------------------|
| Install on Left<br>Part: PTOINDV3585 | Install on Right<br>Part: PTOINDV3586 | Do Not Install<br>Part: PTOINDV3719 |
|--------------------------------------|---------------------------------------|-------------------------------------|

**28.3. Stealth Pelvic/Thigh Guides**

Pelvic/thigh guides may interfere with the joystick receiver block and down posts. Please provide quantity next to check box.

**1. Pelvic/Thigh Guide Bracket**

|                                                                                                          |                                                                                                                              |
|----------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|
| Flip Down Bracket (TWBTLADD-HW) ..... \$236.25 ea.<br>Part: INDPART2804. HCPCS: E1028                    | Flip Down with LEFT Full Surface Contact ..... \$236.25 ea.<br>Feature (TWBADD-L & LPHW407)<br>Part: POS140428. HCPCS: E1028 |
| Removable Bracket (LP430)..... \$262.50 ea.<br>Part: INDPART2805. HCPCS: E1028                           | Flip Down with RIGHT Full Surface Contact..... \$236.25 ea.<br>Feature (TWBADD-R & LPHW407)<br>Part: POS140429. HCPCS: E1028 |
| Removable Bracket with 1" Offset ..... \$400 ea.<br>(LP430 & LHW-121)<br>Part: INDPART2854. HCPCS: E1028 |                                                                                                                              |

**2. Pelvic/Thigh Guide Adapter**

|                                                                                                                                                                                         |                                                                                                                                          |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|
| 2" Height Adapter Plate (LHW-122)..... \$26.25<br>Part: INDPART2821. HCPCS: K0108<br>Height adapter is needed with TRU-Comfort seat cushions. Not available with the flip down bracket. | 2 1/4" Adapter Plate (LHW-129) ..... \$53.50<br>Part: POS140430. HCPCS: K0108<br>Only available with the flip down full surface bracket. |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------|

**3. Pelvic/Thigh Guide Pad**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

|                                                                                        |                                                                                        |
|----------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| Short Pad, 4"x4" (TWBADD-4X4)..... \$75 ea.<br>Part: INDPART2801. HCPCS: E0953^/E0956  | Long Pad, 4"x12" (TWBADD-4X12)..... \$75 ea.<br>Part: INDPART2803. HCPCS: E0953^/E0956 |
| Medium Pad, 4"x8" (TWBADD-4X8)..... \$75 ea.<br>Part: INDPART2802. HCPCS: E0953^/E0956 |                                                                                        |

**4. Pelvic/Thigh Guide Installation**

|                                      |                                       |                                     |
|--------------------------------------|---------------------------------------|-------------------------------------|
| Install on Left<br>Part: PTOINDV3585 | Install on Right<br>Part: PTOINDV3586 | Do Not Install<br>Part: PTOINDV3719 |
|--------------------------------------|---------------------------------------|-------------------------------------|

**28.4. Stealth i-Fit Chest Harness**

|                                                                                       |                                                                                      |
|---------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------|
| Upper Classic, Medium (IFPS1103) ..... \$142 ea.<br>Part: ST-IFPS1103-Q. HCPCS: E0960 | Upper Contour, Medium (IFPS2103)..... \$142 ea.<br>Part: ST-IFPS2103-Q. HCPCS: E0960 |
| Upper Classic, Large (IFPS1104) ..... \$142 ea.<br>Part: ST-IFPS1104-Q. HCPCS: E0960  | Upper Contour, Large (IFPS2104) ..... \$142 ea.<br>Part: ST-IFPS2104-Q. HCPCS: E0960 |

**28.5. Stealth i-Fit Chest Straps**

|                                                                                               |                                                                                                 |
|-----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
| Chest Strap, Medium w/ Belt Anchor (IFCS4113).. \$90 ea.<br>Part: ST-IFCS4113-Q. HCPCS: E0960 | Chest Strap, Large w/ Belt Anchor (IFCS4114)..... \$90 ea.<br>Part: ST-IFCS4114-Q. HCPCS: E0960 |
|-----------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|

**28.6. Stealth Swing-Away Headrest Mounts**


These swing-away headrest mounts are only available with the Stealth Comfort Series headrests. The 1" ball kit must be selected to mount the facial lateral to the headrest bracket. Also select a facial pad(s) with the facial lateral(s).

**1. Facial Lateral Hardware**

|                                                                                                                  |                                                                                                                             |
|------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Facial Lateral for Comfort Series..... \$294.25 ea.<br>Right Side (CPS9-R-HW)<br>Part: INDPART2647. HCPCS: E1028 | Facial Laterals for Comfort Series ..... \$347.75 ea.<br>Left and Right Side (CPS9-B-HW)<br>Part: INDPART2649. HCPCS: E1028 |
| Facial Lateral for Comfort Series..... \$294.25 ea.<br>Left Side (CPS9-L-HW)<br>Part: INDPART2648. HCPCS: E1028  |                                                                                                                             |

**2. Facial Lateral Pads and Accessories**

|                                                                                                                                                        |                                                                                                                                             |
|--------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| 1" Ball Kit for Facial Laterals (TWB-1BK) (D) ..... \$80.25<br>Part: POS150532<br>Required with facial laterals to mount to comfort headrest brackets. | Facial Pad, 4"x2" (900-Q)..... \$53.50 ea.<br>Part: ST-900-Q. HCPCS: K0108                                                                  |
| Egg Switch Mount (ESM) (A) ..... \$52.50 ea.<br>Part: INDPART2815                                                                                      | Facial Pad, 6"x2" (920-Q)..... \$53.50 ea.<br>Part: ST-920-Q. HCPCS: K0108                                                                  |
| Figure 8 Bracket with 13" Rod..... \$267.50 ea.<br>and Egg Switch Mount (B)<br>Part: ST-F8675-13-Q                                                     | Swing-Away Joystick Mount,..... \$250 ea.<br>18" Rod (C)<br>Part: ST-SUS18-Q<br>Allows mounting of facial pads and attaches to the headrest |



Positioning components section continued on following page.

**28. POSITIONING COMPONENTS CONTINUED****28.7. Stealth Elbow Support****1. Elbow Support Bracket**

Left, Swing-Away, Elbow Support Bracket .....\$183.75 ea.  
 Tube Style Arms (ARMS-ES-L)  
 Part: POS124064. HCPCS: E1028

Right, Swing-Away, Elbow Support Bracket.....\$183.75 ea.  
 Tube Style Arms (ARMS-ES-R)  
 Part: POS124066. HCPCS: E1028

**2. Elbow Support Pad**

3"Wx3"H Elbow Pad (ARMS-ES-3X3) .....\$52.50 ea.  
 QTY Part: POS124065. HCPCS: K0108

4"Wx4"H Elbow Pad (ARMS-ES-4X4).....\$52.50 ea.  
 QTY Part: POS124067. HCPCS: K0108

**28.8. Stealth Foot Boxes**

^Please verify the appropriate billing code to be used for this item with all third-party payors.

Stealth Small Foot Boxes, 5"Wx8"Dx7"H.....\$498.75  
 (SPLE107) (pair)  
 Part: INDPART2786. HCPCS: E0954^

Stealth Large Foot Boxes, 7"Wx12"Dx11"H.....\$498.75  
 (SPLE109) (pair)  
 Part: INDPART2788. HCPCS: E0954^

Stealth Medium Foot Boxes, 6"Wx10"Dx9"H ....\$498.75  
 (SPLE108) (pair)  
 Part: INDPART2787. HCPCS: E0954^



Stealth Foot Boxes

**28.9. Stealth Padded Shoe Holder**

Stealth Small Padded Shoe Holder, .....\$266  
 3.75"Wx6.25"D (SPLE111) (pair)  
 Part: INDPART2783. HCPCS: E0951+E0952

Stealth X-Large Padded Shoe Holder, .....\$266  
 5"Wx10.75"D (SPLE114) (pair)  
 Part: INDPART2785. HCPCS: E0951+E0952

Stealth Medium Padded Shoe Holder, .....\$266  
 4"Wx8.25"D (SPLE112) (pair)  
 Part: INDPART2784. HCPCS: E0951+E0952



Stealth Shoe Holders

**29. Q-LOGIC 3 MULTIPLIER HARNESS**

**REQUIRED.** Multiple controllers or power accessories may require an additional multiplier harness or special order. Reference the chart below to determine the appropriate selection. The rows refer to the total amount of arm mounted devices and the columns refer to the total amount of back mounted devices. The online order configurator will automatically make this selection under specialty controls.

No Additional Multiplier (NA)..... No Charge  
 Part: Configuration selection

Secondary Multiplier (SM)..... \$221.21  
 Part: HAR144004

Special Order Configuration (SOF)  
 Part: Configuration selection

|             |                         |
|-------------|-------------------------|
| Arm Mounts  | Joystick (Hand Control) |
|             | Enhanced Display        |
|             | iAccess                 |
|             | USB Charger             |
|             | SAJ                     |
| Back Mounts | Attendant Control       |
|             | SCIM                    |
|             | Sip n Puff              |
|             | ECU                     |

|                             |                         |                            |         |       |       |        |
|-----------------------------|-------------------------|----------------------------|---------|-------|-------|--------|
| NA No additional Multiplier | SM Secondary Multiplier | SOF Special Order Required |         | 1 Arm | 2 Arm | 3+ Arm |
|                             |                         |                            | 0 Back  | NA    | SM    | SOF    |
|                             |                         |                            | 1 Back  | NA    | SM    | SOF    |
|                             |                         |                            | 2 Back  | SM    | SM    | SOF    |
|                             |                         |                            | 3+ Back | SOF   | SOF   | SOF    |

Q-Logic 3 Multiplier with Synergy Seat Chart

Q-Logic 3 Multiplier with Synergy Seat Chart

**30. REFERENCE ORDER FORMS**

**Optional.** Select below to note that additional order forms will accompany this order form. This will also be a link in interactive forms.

[Alternative Controls & Electronics Order Form](#)

[Positioning Components Order Form](#)